

**bbe Moldaenke GmbH**

Preetzer Chaussee 177

24222 Schwentinal

Germany

+49 (0) 431/380 400

+49 (0) 431/380 4010

bbe@mbe-moldaenke.de

**Phone**

**Fax**

**E-mail**

# AlgaeOnlineAnalyser II

## User Manual

Version 2.10 E1, August 2020



---



---

**TABLE OF CONTENTS**


---



---

<b>Table of Contents .....</b>	<b>2</b>
<b>General Safety Precautions .....</b>	<b>9</b>
<b>The AlgaeOnlineAnalyser .....</b>	<b>10</b>
<b>Use of the AlgaeOnlineAnalyser .....</b>	<b>10</b>
<b>Components .....</b>	<b>10</b>
<b>The Sensor Unit with Cleaning Device.....</b>	<b>11</b>
<b>Pump .....</b>	<b>12</b>
Operating .....	12
Loading of the Tube: .....	12
Fitting a Pump Head .....	13
<b>Valves (option) .....</b>	<b>13</b>
<b>Environmental Conditions.....</b>	<b>13</b>
<b>Electrical Connections .....</b>	<b>14</b>
Supply Connection .....	14
Data Interfaces.....	14
RS232 data output .....	14
Ethernet .....	14
Analog Output – 4-20 mA (Option).....	14
Relay Output (Option) .....	15
Connecting the relay output device.....	16
<b>Hydraulic Connections .....</b>	<b>17</b>
Connections of the chlorophyll sensor .....	17
Without integrated sample pump .....	17
With integrated sample pump.....	17
Multichannel version .....	18
Tubing diagram with valves and pump.....	18
Tubing diagram with valves only .....	19
Drain .....	19
<b>Quick Start .....</b>	<b>19</b>
<b>Measurement Procedures .....</b>	<b>19</b>
Measuring Sequence (Standard sequence) .....	19
Measuring Sequence (flow through) .....	19
Measuring Sequence (multichannel version) .....	20
Interval pumping vs. continuous pumping .....	20
<b>bbe++ Software.....</b>	<b>21</b>

<b>Software Set-Up .....</b>	<b>21</b>
Desktop icon of the bbe++ software .....	23
bbe++ software in the Start menu .....	23
<b>Software Update .....</b>	<b>24</b>
<b>Structure of the software .....</b>	<b>24</b>
<b>Menu structure (“Basic functionality”) .....</b>	<b>25</b>
<b>Menu structure (“Advanced functionality”) .....</b>	<b>26</b>
<b>General considerations for storing and displaying data and parameters in bbe++ .....</b>	<b>27</b>
What is stored in a bbe++ database? .....	27
Where is the database stored? .....	27
How to display the data? .....	27
Simultaneous display of data of different types of instruments .....	27
Who can access to the parameters of the instruments? .....	28
Desktop of the bbe++ Software. ....	28
Toolbar .....	29
Toolbar instrument .....	29
Graphic- und table views .....	29
Serial number, voltage of the battery, internal clock .....	29
<b>The Menus .....</b>	<b>30</b>
File Menu .....	30
New .....	30
Open .....	31
Close .....	31
Import (FluoroProbe / AlgaeTorch / BenthosTorch) .....	31
Export (ASCII) .....	33
Export (Database) .....	33
Export (KMZ) – AlgaeTorch / BenthosTorch / FluoroProbe .....	33
Export (ToxTest protocol) – AlgaeLabAnalyser only .....	33
Summary View .....	33
1. ... 2. ... 3. ....	34
Data Directory .....	35
Exit .....	35
Edit Menu .....	35
Copy .....	35
Edit Comment – advanced functionality only .....	36
Delete Data – advanced functionality only .....	37
Parameters Menu .....	38
Overview .....	40

---

History – advanced functionality only .....	41
Templates .....	44
Online.....	46
Recalc – advanced functionality only.....	48
Matrix view .....	51
View Menu.....	52
Toolbar .....	52
Status Bar.....	52
Measurement .....	52
Macro.....	53
Data Navigation .....	53
Period.....	53
Session .....	53
Connection State .....	53
Auto-Update .....	54
Application look .....	54
Window – Advanced functionality only .....	54
New Table Window.....	54
New XYGraph Window .....	54
New Single Value View .....	55
Cascade .....	56
Tile Horizontally/Vertically .....	57
Select Data .....	57
View Editor.....	58
Comment Fields .....	73
1. ... 2. ... 3. ....	74
Tools.....	74
Settings .....	74
Enable advanced/basic functionality .....	83
Connect COM port .....	83
Device .....	84
AOA device Submenu.....	84
SQL Script.....	88
Calibration .....	88
Sample temperature.....	88
Transmission .....	88
Chlorophyll-A .....	89
Offsets.....	89

---

Calibration solution for the „ultra filtrated water offset and the yellow substances.....	90
Offsets (ultrafiltrated).....	90
Offsets (distilled).....	91
Fingerprints.....	92
Correction factors.....	93
Correction of a fingerprint.....	93
Calibration of Fingerprints with Subsequent Wet-Chemical Analysis:.....	93
Help.....	94
About.....	94
bbe++ manuals.....	94
<b>The Views.....</b>	<b>95</b>
Table views.....	95
History of Parameters.....	96
Edit comment.....	96
Delete data.....	97
Data Fields Selection.....	98
Sorting of the list view.....	98
XY graph views.....	99
Select data of interest.....	100
Right-click menu.....	102
Scaling of the axis.....	102
Autoscaling.....	104
Parameters.....	104
Default scaling.....	104
<b>Parameters.....</b>	<b>105</b>
Process parameters.....	105
LED measuring time.....	105
Measurement interval.....	105
Number of chlorophyll (f) measurements.....	105
Measurement duration (f, fm, fo, transmission).....	105
Minimum chlorophyll (f) concentration.....	105
Enter comment at start of measurement.....	106
Activity.....	106
Continuous measurements.....	106
Number of available sample channels.....	106
Pump mode.....	106
Pump time start up.....	106
Adaptation time (f, fo).....	107

---

Integration time fm measurement .....	107
Alarm parameters .....	107
Alarm limit .....	108
Enable alarm .....	108
Common parameters .....	108
Warm-up time .....	109
Serial number .....	109
Software version .....	109
Fit parameters .....	109
Name of class .....	109
Enabled for fit .....	109
Yellow Substances Correction .....	110
Offsets .....	110
Transmission offsets / gradients .....	111
Fingerprints .....	111
Standard deviations .....	112
Date of calibration .....	112
Global correction factor .....	112
Cell factors .....	112
Fo measurement brightness .....	112
Measurement parameters .....	113
DA values of the LEDs .....	113
Temperature and Humidity offsets / gradients .....	113
Required values of the LEDs .....	113
Temperature correction .....	114
<b>Operating Principles .....</b>	<b>115</b>
<b>Principles of Operation .....</b>	<b>115</b>
Determination of different algae .....	115
Determination of yellow substances .....	115
Data analysis .....	115
Offset .....	115
Algae classes .....	116
Measurement of the Algae Activity .....	116
Measurement of Further Values and their Calculation .....	117
In General .....	117
Transmission .....	117
Temperature of the Sensor .....	118
Brightness of the LED .....	118

---

<b>Measuring Procedure</b> .....	<b>118</b>
Initial Measuring Procedure .....	118
Offsets.....	118
<b>Measuring Procedure for Low Chlorophyll Concentration</b> .....	<b>118</b>
Background Fluorescence .....	119
Cleaning .....	119
Measuring Time .....	119
Sedimentation .....	119
Sample Preparation .....	119
<b>Calibration Procedure</b> .....	<b>119</b>
General Considerations: .....	119
Calibration of the Offsets:.....	119
Calibration of the Algae Fingerprints:.....	119
Calibration of the Yellow Substances: .....	120
Calibration by Factors: .....	120
Calibration of Fingerprints with Subsequent Wet-Chemical Analysis: .....	120
<b>Batch Mode Operation</b> .....	<b>120</b>
<b>Data Evaluation</b> .....	<b>122</b>
<b>Serial data exchange</b> .....	<b>123</b>
<b>Technical Data</b> .....	<b>125</b>
<b>Maintenance</b> .....	<b>127</b>
<b>Weekly Maintenance</b> :	<b>127</b>
Maintenance Work that Has to be Done: .....	127
<b>Monthly Maintenance</b> .....	<b>127</b>
Maintenance Work that Has to be Done: .....	127
<b>Hints for Maintenance Work</b> .....	<b>127</b>
Cleaning of the Hose System .....	127
Exchange the Pump Hoses.....	127
Hose diameter and hose types .....	127
Clean the Measuring Cuvette .....	128
<b>Change Fuses</b> .....	<b>129</b>
<b>Fault finding</b> .....	<b>130</b>
<b>The Algae Online Analyser will not Start</b> .....	<b>130</b>
<b>Pump is not running</b> .....	<b>130</b>
<b>No Algae Concentration Shown</b> .....	<b>130</b>
<b>Genty Results are Zero</b> .....	<b>130</b>
<b>Results Too Low</b> .....	<b>130</b>

**Results Too Low or Too High ..... 130**

**Humidity Detected in Sensor ..... 130**

***Troubleshooting* ..... 131**



## GENERAL SAFETY PRECAUTIONS

---

---

The bbe AlgaeOnlineAnalyser is an instrument for the determination of the concentration of chlorophyll in water. It should only be used for this purpose.

- only instructed staff should operate this instrument.
- consult appropriate safety manual in case of using hazardous compounds and solutions. Wear gloves, coat and safety goggles.
- electrical connection of the instrument should only be carried out by trained staff.
- the instrument should only be opened by authorized staff.
- please refer to the safety instructions of any chemicals used with the AlgaeOnlineAnalyser.
- changes to electrical connections and circuits may cause damage to the instrument and lead to loss of warranty.

---

## THE ALGAEONLINEANALYSER

---

### USE OF THE ALGAEONLINEANALYSER

---

The AlgaeOnlineAnalyser is an instrument for the analysis of algae in rivers, lakes or reservoirs. It is usually used in a monitoring station. Depending on the configuration of the instrument, there are various parameters which can be determined:

- the total concentration of chlorophyll
- the concentration of up to 5 algae groups (algae differentiation)
- the activity or production rate of these algae groups (Option)
- the transmission of 5 wavelengths (Option)
- detection of yellow substances


The user receives an instrument for the complete survey of the algae content of waterways.

### COMPONENTS

---

The main components of the bbe AlgaeOnlineAnalyser are:

- The sensor unit with cleaning device
- Housing with power-supply and pump (optional valves for multichannel operation)
- PC unit with touchscreen display




AlgaeOnlineAnalyser front view

## THE SENSOR UNIT WITH CLEANING DEVICE

Here the actual measurement is processed. The sample is pumped into the cuvette from below. There the fluorescence is measured with light of different intensity and colour. The data is analysed by a microcontroller integrated into the sensor unit and transferred to the PC. After a certain number of measurements the cuvette is cleaned automatically by the cleaning device.

The humidity detector prevents greater damage to the sensor in case of leakage.



cross-section of the sensor unit (Chlorophyllsensor)

PUMP

The bbe Algae Online Analyser also uses a peristaltic pump – the only degradable part is the hose.


The rolls on the rotor press the hose against the hose saddle and close it at that point. So there is a constant volume between the two rolls that is delivered in the direction of the rotor.

The quantity delivered depends on the number of revolutions of the pump and on the inner diameter of the hose.

Another advantage is the use of several channels per head. This offers the possibility to have the pumping rates depend on each other and to be determined by the relation of the inner diameters of the hoses.


The hoses of the sample pump are put under considerable strain so we recommended the hoses be changed every 4 weeks.

OPERATING





- 1 Tube clamp setting.** Predisposizione del morsetto del tubo. Justera slangklämmen. Einstellung der Schlauchklemme. Fjar le abrazadera para tubos. Afstellen van slangklemmen. Règlez de serrage du tube. Predisposizione del morsetto de tubo. Putken kiinnikkeen asettaminen. Innstilling av slang klemme. Innstilling af slangklemmer.
- 2 Tube wall thickness.** Spessore delle pareti del tubo. Slangväggens grovlek. Wandstärke des Schlauchs. Grosor de la pared del tubo. Wanddicke pomslang. Epaisseur de la paroi du tube. Spessore delle pareti del tubo. Putken seinämän vahvuus. Vegtykkelse på slange. Vægtykkelse.
- 3 Set for tube bore size.** Predisposto per le dimensioni del diametro del tubo. Ställ in för olika slangstorlekar. Satz für Größe der Schlauchbohrung. Fjar para el tamaño di-metro interior del tubo. Insteloid voor slangdoorliten. Règler le diamètre interne du tube. Predisposto per le dimensioni del diametro del tubo. Asetus putken halkaisijan koolle. Innstilt størrelse på slange diameter. Innstilling af slangediameter

Opening the pumphead. Apertura della testa della pompa. Öppna pumphuvudet. Öffnen des Pumpenkopfs. Abrir el cabezal. Openen van pompkop. Ouverture de la tête de pompe. Apertura della testa della pompa. Pumpun pään avaaminen. Åpning av pumpehode. Åbning af pumpehovedet.



LOADING OF THE TUBE:



**FITTING A PUMP HEAD****VALVES (OPTION)**

In case of a multichannel instrument, pinch valves are used. Please use the original type of tubes for proper operation.

**ENVIRONMENTAL CONDITIONS**

For the installation of the bbe Algae Online Analyser a frost-free room has to be chosen. The temperature must not exceed 35°C.

Parameters	Values
Temperature range during use	5 - 35°C
Temperature range during storage	0 - 50°C
Relative humidity	up to 95 %, non-condensed
Exposure to sunlight	No direct exposure to sunlight
Minimum distance from the wall	10 cm
Maximum pressure	0.5 bar above atmospheric pressure

A free drain should also be provided. The bbe Algae Online Analyser is built as a tabletop instrument.

## ELECTRICAL CONNECTIONS

---

### SUPPLY CONNECTION

The bbe AlgaeOnline Analyser is supplied with a mains cable with a safety plug. When connecting the instrument, please ensure a waterproof connecting socket.

### DATA INTERFACES

#### RS232 data output

##### **Important:**

**Please switch off the AlgaeOnlineAnalyser and the external device before connecting the RS232 cable.**

The AOA provides a RS232 data output. This is available on a free RS232 port of the external PC or via the RS232 output cable when using an internal PC.

#### Ethernet


The AOA can be integrated in a LAN. The LAN cable has to be connected directly to the PC. To do so the cable has to be lead through the cable inlet at the right side of the housing.



#### ANALOG OUTPUT – 4-20 mA (OPTION)


Algae Online Analysers can be equipped with 2 or more analog outputs as an option. In this case, there are one or more 4-pin plug connector(s) mounted on the instrument. An extra plug that can be attached to the on-site connection cable is also supplied.

Please note: Common ground lines are not supported. The Ground lines of the 4-20mA outputs must be isolated from each other.



1: output 0 +

2: output 0 -




connector  
from the  
'screw-side'

Gnd: output 1 -

3: output 1 +

1: output 2 +

2: output 2 -



connector  
from the  
'screw-side'

Gnd: output 3 -

3: output 3 +

In case of more analog outputs, the following connectors are marked consecutively.

### RELAY OUTPUT (OPTION)

The Algae Online Analyser can be equipped with a relay output device.

### Connecting the relay output device

There are two (or more) 7-pole connectors on the right side of the instrument. The corresponding cable connector is delivered together with the accessories.




Table of connection:

Connector	Pole <i>Printed on the connector</i>	Function
1	1	Relay 2 / Common
1	2	--
1	3	Relay 2 / Normally Open
1	4	Relay 1 / Common
1	5	--
1	6	Relay 1 / Normally Open
2	1	Relay 4 / Common
2	2	Relay 4 / Normally Closed
2	3	Relay 4 / Normally Open
2	4	Relay 3 Common
2	5	Relay 3 Normally Closed
2	6	Relay 3 Normally Open




---

## HYDRAULIC CONNECTIONS

---

### CONNECTIONS OF THE CHLOROPHYLL SENSOR



#### WITHOUT INTEGRATED SAMPLE PUMP

The sample has to flow continuously through the instrument. For this, a maximum pressure of 0.2 bar is allowed. Please ensure that this pressure is not exceeded, even if the effluent is blocked. The inlet is on the bottom in the center of the removable lid.

#### WITH INTEGRATED SAMPLE PUMP


Connect the intake of the integrated pump to the water source and the outlet of the pump to the connector at the bottom of the sensor of the AlgaeOnlineAnalyser.


### MULTICHANNEL VERSION

Integrating valves in the AOA extends the instrument to a multi-channel device. The valve can be used with and without the peristaltic pump. The pump is required as soon as the pressure of one of the sample streams is not sufficient.



### Tubing diagram with valves and pump



Tubing diagram with valves onlyDRAIN

The drain of the instrument needs to be fed into an open, pressure-free drain.

QUICK START

- Connect the mains
- Connect the LAN / RS232 / 4-20mA
- Connect the hoses
- Close the peristaltic pump
- Switch the instrument on and let the PC boot
- Adjust pumping time and procedure according to your needs
- Start the measurement using the “start” button
- The measurement can be stopped at any time via the “stop” button

MEASUREMENT PROCEDURES

First of all, a survey of the single steps of a measurement.

MEASURING SEQUENCE (STANDARD SEQUENCE)

1. Pumping during startup or opening the sample valve
2. Brightness tuning of the internal LED
3. Pumping sample or opening the sample valve
4. Measurement (including activity measurement if available)
5. Cleaning
6. Waiting for the next measurement
7. Go to step 2

MEASURING SEQUENCE (FLOW THROUGH)

1. Flushing during startup
2. Brightness tuning of the internal LED
3. Measurement (concentration measurement only)
4. Cleaning
5. Waiting for the next measurement
6. Go to step 2

### MEASURING SEQUENCE (MULTICHANNEL VERSION)

1. Pumping channel 1 during startup or opening the sample valve 1
2. Pumping channel 2 to n during startup or opening the sample valve 2 to n
3. Brightness tuning of the internal LED
4. Pumping sample or opening the sample valve for the current channel
5. Measurement (including activity measurement if available)
6. Cleaning
7. Go to step 2 for channel 1 to n
8. Waiting for the next measurement

### INTERVAL PUMPING VS. CONTINUOUS PUMPING

Two operation modes for the pump and/or the valves are selectable:

- Interval mode or
- Continuous mode.

#### Interval mode

The pump is working (or the valve is opening) before each measurement and it is not working during the measurement.

Advantages;

This enables to perform activity measurements.

The working time of the pump and the tube is short, the lifetime extended.

Disadvantages;

Measurement is slower

In multichannel applications it takes longer to exchange the sample

#### Continuous mode.

The pump is working (or the valve is opening) all the time for fast response. Not suitable for activity measurements.

---

## BBE++ SOFTWARE

---

The bbe++ software is delivered together with bbe instruments. If a PC is delivered with the instrument or integrated in the instrument, the software is already installed.

It provides the following functions:

- operation, control and calibration of bbe instruments
- data analysis and display in tables and diagrams
- export in different formats

This chapter describes the general functions of the bbe++ software. The examples may show data or parameters of other bbe instruments. The data and parameters of your instrument can easily be handled accordingly. Whenever this manual deals with special features of a certain instrument, this will be explained in the special instrument chapter.

The data and parameters of the instruments are stored in a database. A single database may contain the data of different instruments of the same type and also different types of instruments.


## SOFTWARE SET-UP

---


The software is stored on a CD delivered with the instrument. The latest version can also be downloaded from the bbe website after registration (<http://www.bbe-moldaenke.de/log-in/>).

The Windows autorun function should start the installation automatically. If it does not, open the Explorer and click on **setup.exe** on the CD for installation.


Welcome window: click "NEXT >"




Choose the destination folder and click "NEXT >"




Choose a new start menu folder and click "NEXT >"



Choose your preferred language and click "NEXT >"



Final success window of the bbe++ installation. Click "FINISH".



### DESKTOP ICON OF THE BBE++ SOFTWARE


The setup installs an icon for the bbe++ software on the desktop (if chosen during installation).



### BBE++ SOFTWARE IN THE START MENU

The bbe++ folder in the start menu contains:

- bbe++ software
- installed manuals
- link to the data folder of bbe++



## SOFTWARE UPDATE

---

Updates of the bbe++ software and manuals can be downloaded from the bbe website <http://www.bbe-moldaenke.de> after registration.

## STRUCTURE OF THE SOFTWARE


---

The following menus are available in the bbe++ software. Whether a menu is displayed, depends on the access level set.


- FILE contains all the input/output functions.
- EDIT contains functions to copy data and graphics.
- PARAMETER contains functions to set the parameters of an instrument and to adapt the parameters of data already measured.
- VIEW contains functions to show and hide toolbars.
- WINDOW contains functions to show the data in different types of tables and graphics as well as editors to change the layout (advanced functionality level only)
- TOOLS contains options to change settings that influence the behaviour of the bbe++ software and the way of operating the instruments.
- CALIBRATION contains the items to calibrate different bbe instruments (advanced functionality level only)
- HELP contains information about the current bbe++ version.



## MENU STRUCTURE (“BASIC FUNCTIONALITY”)



MENU STRUCTURE ("ADVANCED FUNCTIONALITY")



---

## GENERAL CONSIDERATIONS FOR STORING AND DISPLAYING DATA AND PARAMETERS IN BBE++

---

### WHAT IS STORED IN A BBE++ DATABASE?


The measuring results and instrument parameters of different bbe instruments are stored. It is possible to retrieve the corresponding instrument parameters for each single result.

In the so-called "Summary View" the serial numbers of the instruments and the time ranges for all data stored in this database are shown. This view can be used to select the desired data as well as the corresponding parameters for display.

The database files are marked with the extension \*.bdb (bbe database)

### WHERE IS THE DATABASE STORED?

The database of bbe++ is stored in the standard application folder of Windows. The name of the folder depends on the operation system and the language. For example: C:\ProgramData\bbe++ in the Windows 7 English version. To access this folder easily, there is a link in the file menu of bbe++ or start menu of Windows called "Data directory":



### HOW TO DISPLAY THE DATA?

In the bbe++ software, different views have been pre-defined for each instrument:

- graphics
- tables
- single dataset
- data export

User-defined views may be added as well.


### SIMULTANEOUS DISPLAY OF DATA OF DIFFERENT TYPES OF INSTRUMENTS

After opening a database file, the data of different instruments of the same type may be displayed simultaneously. To display the data of another type of instrument from the same database file, this file can be opened more than once.

WHO CAN ACCESS TO THE PARAMETERS OF THE INSTRUMENTS?

The bbe++ software has 2 access levels – basic functionality and advanced functionality. Depending on the different access levels, more or fewer menu items and parameters are shown. There are different settings for reading or writing access. The advanced functionality level can be password protected.

DESKTOP OF THE BBE++ SOFTWARE.



Toolbar

Toolbar Instrument

Data table

Graph 1

Graph 2

Serial number, voltage of the battery, internal clock

### TOOLBAR

Toolbar to navigate within the datasets.

### TOOLBAR INSTRUMENT

Toolbar to operate the instrument

### GRAPHIC- UND TABLE VIEWS

User configurable windows with graphs and tables.


### SERIAL NUMBER, VOLTAGE OF THE BATTERY, INTERNAL CLOCK

Display of the data of the connected instrument: serial number, voltage of the battery and internal clock.

## THE MENUS


### FILE MENU

The File Menu contains all the input/output functions.




### New


Creates a new database file for data and parameters.



If the new database already exists, the existing database is deleted after confirming the following message:




The new and empty database is displayed.



New data can be entered in the database by importing or by starting a measurement.

### Open

Opens an existing database. If there are already data in the database, the "Summary View" window appears. The type and serial number of the instrument can be chosen as well as the time range. The following example shows a database containing data from different FluoroProbes.




### Close

Closes the active database.

### Import (FluoroProbe / AlgaeTorch / BenthosTorch)


This function is used to import \*.FLP files into the database. These files can be results transferred from the FluoroProbe to a USB stick or files generated by the first version of the FluoroProbe software.

To import, please choose the file from the file system:



It is possible to select more than one file by pressing the shift key while selecting the file.

FluoroProbe only: Afterwards choose how to apply the air pressure and hence the calculation of depth.



#### Air pressure as given in file (FluoroProbe and AlgaeTorch 100 only)

Air pressure was measured at the beginning of the measurement. This value is used.


#### Air pressure from first data set (FluoroProbe and AlgaeTorch 100 only)

The first measurement of air pressure took place in air. This is the appropriate selection if an autostart-plug without a PC was used (for further information see the FluoroProbe-manual).

#### Air pressure (FluoroProbe and AlgaeTorch 100 only)


Allows the entry of a constant value. To be used if the FluoroProbe was submerged for a longer time and unable to measure the air pressure by itself.

In case of FluoroProbe data, the results are recalculated from the raw data after importing.



The calculation procedure takes place in subsequent steps. After importing all data, a success message is shown.





### Export (ASCII)

This item exports the data as an ASCII file to load them into other programs. To select columns to be exported and the separator that is used, please see the description of the View Editor. The view editor allows the user to define different export formats for each instrument. All export formats of the current instrument are shown in the menu:



The default export view is “Export all”.

After clicking on “Export all”, the name and folder of the exported file can be selected.

Hint: To transfer data to Excel, use the copy and paste feature in the “Edit” menu.

### Export (Database)

This item is used to export parts of a database to a new (and smaller) one. The exported database can be read with bbe++ again. This might be useful when parts of the database are transferred via email.

The currently selected data (“WINDOW → SELECT DATA”) are exported.

### Export (KMZ) – AlgaeTorch / BenthosTorch / FluoroProbe


This item is used to export geo data from instruments with GPS receiver to be displayed in Google Earth.

### Export (ToxTest protocol) – AlgaeLabAnalyser only

This item is used to export the results of toxicity tests conducted with the AlgaeLabAnalyser to Excel.

### Summary View

Opens a new window that gives an overview of all the data and instruments in the current database:




The database bbe++.bdb above contains data of two types of instruments: AlgaeTorch and FluoroProbe.

The serial numbers and the point in time of the first measurement and the last measurement stored in the database are shown below the heading with the name of the instrument.

The header of the window (bbe++.bdb [FP-18-03] :3 has the following meaning:

- the name of the database bbe++.bdb
- the serial number of the instrument  
for which the data in the data windows are currently shown [FP-18-03]
- the number of the window if there is more than one open window :3

The summary view can be used to directly access the data or parameters of one of the listed instruments. To do so, use the right mouse click on the serial number:



A selection of different periods and parameters are shown. After selecting one of the given periods, it is shown using the default view.

The screenshot shows a window titled "bbe++.bdb [FP-15-02]". It displays a table with the following data:

Date/Time [date]	Conc 0 [Green Algae] [µg/l]	Conc 1 [Bluegreen] [µg/l]	Conc 2 [Diatoms] [µg/l]
22.02.2007 15:27:25	0	1,34	
22.02.2007 15:27:34	0	1,43	
22.02.2007 15:27:40	0	1,41	
22.02.2007 15:27:47	0	1,37	
22.02.2007 15:27:53	0	1,45	
22.02.2007 15:28:00	0	1,4	
22.02.2007 15:28:06	0	1,34	
22.02.2007 15:28:13	0	1,35	
22.02.2007 15:28:20	0	1,38	
22.02.2007 15:28:26	0	1,42	
22.02.2007 15:28:35	0	1,3	
22.02.2007 15:28:41	0	1,47	
22.02.2007 15:28:48	0	1,4	
22.02.2007 15:28:54	0	1,33	
22.02.2007 15:29:01	0	1,46	
22.02.2007 15:29:08	0	1,39	
22.02.2007 15:29:14	0	1,37	

For further information about the parameters window, please see the chapter Parameters.

1. ... 2. ... 3. ...

Names of the last database files opened for quick reload.

## Data Directory


Opens the explorer with the current data directory to copy or rename the data base. This is the same function as in the start menu up to Windows 7.

## Exit

Terminates the program.

## EDIT MENU

This menu contains commands to copy data/graphics from bbe++ and add them to the clipboard.



## Copy

The content of the active window will be copied to the clipboard. Depending on the type of the active window, this is either a graph or a table.

Tables can be pasted to EXCEL directly.

## How to copy data to EXCEL:

Use "WINDOW → SELECT DATA" to select the data you want to copy to EXCEL. Use the list window type that contains all the columns to be exported. Please note: only the columns from the active view are exported.

Date/Time [date]	Conc 0 [Green Algae] [µg/l]	Conc 1 [Bluegreen] [µg/l]	Conc 2 [Diatoms] [µg/l]	Conc 3 [Cryptophyta] [µg/l]
15.06.2009 09:10:50	0	0,813	0	0,495
15.06.2009 09:10:53	0	1,04	0	0
15.06.2009 09:10:55	0	0	0	0
15.06.2009 09:10:57	0	2,28	0	0
15.06.2009 09:11:00	1,08	0,829	0	0
15.06.2009 09:11:02	0	0,142	0	2,09
15.06.2009 09:11:04	0,0475	0,956	0	0,908
15.06.2009 09:11:07	0	0,932	0	0
15.06.2009 09:11:09	0	1,2	0	0
15.06.2009 09:11:12	0	0,74	0	0
15.06.2009 09:11:14	0	0,561	0	0,936
15.06.2009 09:11:16	0	0,698	0	0
15.06.2009 09:11:19	0,945	0,779	0	0
15.06.2009 09:11:21	0	0,804	0	0
15.06.2009 09:11:23	3,35	0	0	0
15.06.2009 09:11:26	1,07	0,488	0	0
15.06.2009 09:11:28	0	0,261	0	2,03
15.06.2009 09:11:31	1,52	0,227	0	0
15.06.2009 09:11:33	0	1,1	0	0
15.06.2009 09:11:35	0	1,58	0	1,97
15.06.2009 09:11:38	5,41	0	0	0,883
15.06.2009 09:11:40	3,26	0	2,18	0
15.06.2009 09:11:42	0,792	0	2,82	0
15.06.2009 09:11:45	2,82	0	1,22	0,85

- click on "EDIT → COPY"
- open Excel with an empty datasheet
- click on "EDIT → PASTE" in Excel

	A	B	C	D	E	F
1	Date/Time	Conc 0 [Green Algae]	Conc 1 [Bluegreen]	Conc 2 [Diatoms]	Conc 3 [Cryptophyta]	Conc 4 [Other]
2	date	µg/l	µg/l	µg/l	µg/l	µg/l
3	15.06.2009 09:10	0	0,813	0	0,495	0
4	15.06.2009 09:10	0	1,04	0	0	0
5	15.06.2009 09:10	0	0	0	0	0
6	15.06.2009 09:10	0	2,28	0	0	0
7	15.06.2009 09:11	1,08	0,829	0	0	0
8	15.06.2009 09:11	0	0,142	0	2,09	0
9	15.06.2009 09:11	0,0475	0,956	0	0,908	0
10	15.06.2009 09:11	0	0,932	0	0	0
11	15.06.2009 09:11	0	1,2	0	0	0
12	15.06.2009 09:11	0	0,74	0	0	0
13	15.06.2009 09:11	0	0,561	0	0,936	0
14	15.06.2009 09:11	0	0,698	0	0	0
15	15.06.2009 09:11	0,945	0,779	0	0	0
16	15.06.2009 09:11	0	0,804	0	0	0
17	15.06.2009 09:11	3,35	0	0	0	0
18	15.06.2009 09:11	1,07	0,488	0	0	0
19	15.06.2009 09:11	0	0,261	0	2,03	0
20	15.06.2009 09:11	1,52	0,227	0	0	0
21	15.06.2009 09:11	0	1,1	0	0	0
22	15.06.2009 09:11	0	1,58	0	1,97	0
23	15.06.2009 09:11	5,41	0	0	0,883	0
24	15.06.2009 09:11	3,26	0	2,18	0	0
25	15.06.2009 09:11	0,792	0	2,82	0	0
26	15.06.2009 09:11	2,92	0	1,33	2,85	0
27	15.06.2009 09:11	6,68	0	0	1,54	0
28	15.06.2009 09:11	6,56	0	0,939	0,492	0
29	15.06.2009 09:11	0	0,111	1,16	0	0
30	15.06.2009 09:11	0	0,363	0	0,598	0
31	15.06.2009 09:11	4,4	0	0	1,77	0

Now the data can be used for all kinds of calculations within Excel.

### Edit Comment – advanced functionality only

To edit the comment of one or more datasets mark the datasets in the table view. A dataset can be marked with a mouse click. Multiple datasets can be marked by pressing Ctrl key while clicking on the dataset required.

Date/Time [date]	Comment	Green Algae [µg/l]	Bluegreen [µg/l]	Diatoms [µg/l]	Cryptophyta [µg/l]	#5 [µg/l]	Yellow substances [µg/l]	Total conc. [µg/l]
27.12.2011 08:41:23		0,33	1,87	0,00	0,00	0,00	0,00	2,20
27.12.2011 08:41:24		0,06	1,58	0,51	0,00	0,00	0,00	2,14
27.12.2011 08:41:24		0,07	1,59	0,00	0,00	0,00	0,00	1,65
27.12.2011 08:41:26		0,00	1,66	0,00	0,00	0,00	0,00	1,66
27.12.2011 08:41:29		0,00	2,37	0,00	0,00	0,00	0,20	2,37
27.12.2011 08:41:30		0,91	1,69	0,02	0,00	0,00	0,00	2,62
27.12.2011 08:41:34		0,00	0,00	0,38	0,82	0,00	0,00	1,20
27.12.2011 08:41:36		0,00	0,25	0,00	0,00	0,00	1,47	0,25
27.12.2011 08:41:39		0,00	0,62	0,00	1,29	0,00	0,00	1,91
27.12.2011 08:41:41		0,00	3,86	0,00	1,24	0,00	0,57	5,10
27.12.2011 08:41:50		0,00	0,00	0,00	2,26	0,00	0,00	2,26
27.12.2011 08:41:51		0,00	0,00	0,00	4,87	0,00	0,00	4,87
27.12.2011 08:41:57		0,11	0,00	2,64	0,46	0,00	0,00	3,21
27.12.2011 08:41:58		0,24	0,00	2,67	0,13	0,00	0,27	3,04
27.12.2011 08:41:59		1,34	0,00	2,64	0,00	0,00	0,00	3,98
27.12.2011 08:41:59		0,00	0,68	3,78	0,00	0,00	0,09	4,46
27.12.2011 08:41:59		1,23	0,00	1,60	1,23	0,00	0,00	4,05
27.12.2011 08:42:00		2,90	0,00	1,36	0,00	0,00	0,09	4,26
27.12.2011 08:42:00		0,95	0,00	2,92	0,14	0,00	0,14	4,01
27.12.2011 08:42:01		0,38	0,00	3,08	0,00	0,00	0,51	3,46
27.12.2011 08:42:01		0,48	0,00	3,57	0,00	0,00	0,00	4,05
27.12.2011 08:42:01		1,63	0,00	2,79	0,00	0,00	0,00	4,43
27.12.2011 08:42:02		0,71	0,00	3,36	0,36	0,00	0,00	4,43
27.12.2011 08:42:02		1,51	0,00	3,02	0,00	0,00	0,00	4,54
27.12.2011 08:42:03		1,53	0,00	2,92	0,00	0,00	0,21	4,45

Go to Edit → Edit comments or use right click on the marked datasets:

	0,07	1,59
	0,00	1,66
	0,00	2,37
History of Parameters		1,69
Edit Comment		0,00
Delete Data		0,25
Data Fields Selection		0,62
	0,00	3,86
	0,00	0,00
	0,00	0,00
	0,11	0,00

Enter the new comment for the selected data:

Date/Time [date]	Comment	Green Algae [µg/l]	Bluegreen [µg/l]	Diatoms [µg/l]	Cryptophyta [µg/l]	#5 [µg/l]	Yellow substances [µg/l]	Total conc. [µg/l]	Transmissior
27.12.2011 08:41:23		0,33	1,87	0,00	0,00	0,00	0,00	2,20	9
27.12.2011 08:41:24		0,06	1,58	0,51	0,00	0,00	0,00	2,14	9
27.12.2011 08:41:24		0,07	1,59	0,00	0,00	0,00	0,00	1,65	9
27.12.2011 08:41:26		0,00	1,66	0,00	0,00	0,00	0,00	1,66	10
27.12.2011 08:41:29		0,00	2,37	0,00	0,00	0,00	0,20	2,37	9
27.12.2011 08:41:30		0,91	1,69	0,02	0,00	0,00	0,00	2,62	9
27.12.2011 08:41:34		0,00	0,00	0,38	0,82	0,00	0,00	1,20	9
27.12.2011 08:41:36		0,00	0,25	0,00	0,00	0,00	1,47	0,25	10
27.12.2011 08:41:39									10
27.12.2011 08:41:41									9
27.12.2011 08:41:50									10
27.12.2011 08:41:51									9
27.12.2011 08:41:57									10
27.12.2011 08:41:58									10
27.12.2011 08:41:59									10
27.12.2011 08:41:59									10
27.12.2011 08:41:59									10
27.12.2011 08:42:00		2,90	0,00	1,36	0,00	0,00	0,09	4,26	10
27.12.2011 08:42:00		0,95	0,00	2,92	0,14	0,00	0,14	4,01	10
27.12.2011 08:42:01		0,38	0,00	3,08	0,00	0,00	0,51	3,46	10
27.12.2011 08:42:01		0,48	0,00	3,57	0,00	0,00	0,00	4,05	10
27.12.2011 08:42:01		1,63	0,00	2,79	0,00	0,00	0,00	4,43	10
27.12.2011 08:42:02		0,71	0,00	3,36	0,36	0,00	0,00	4,43	10
27.12.2011 08:42:02		1,51	0,00	3,02	0,00	0,00	0,00	4,54	10
27.12.2011 08:42:03		1,53	0,00	2,92	0,00	0,00	0,21	4,45	10

Click OK:

Date/Time [date]	Comment	Green Algae [µg/l]	Bluegreen [µg/l]	Diatoms [µg/l]	Cryptophyta [µg/l]	#5 [µg/l]	Yellow substances [µg/l]	Total conc. [µg/l]	Transmissior
27.12.2011 08:41:23		0,33	1,87	0,00	0,00	0,00	0,00	2,20	9
27.12.2011 08:41:24		0,06	1,58	0,51	0,00	0,00	0,00	2,14	9
27.12.2011 08:41:24		0,07	1,59	0,00	0,00	0,00	0,00	1,65	9
27.12.2011 08:41:26	This is the new comm...	0,00	1,66	0,00	0,00	0,00	0,00	1,66	10
27.12.2011 08:41:29	This is the new comm...	0,00	2,37	0,00	0,00	0,00	0,20	2,37	9
27.12.2011 08:41:30	This is the new comm...	0,91	1,69	0,02	0,00	0,00	0,00	2,62	9
27.12.2011 08:41:34	This is the new comm...	0,00	0,00	0,38	0,82	0,00	0,00	1,20	9
27.12.2011 08:41:36	This is the new comm...	0,00	0,25	0,00	0,00	0,00	1,47	0,25	10
27.12.2011 08:41:39	This is the new comm...	0,00	0,62	0,00	1,29	0,00	0,00	1,91	9
27.12.2011 08:41:41	This is the new comm...	0,00	3,86	0,00	1,24	0,00	0,57	5,10	9
27.12.2011 08:41:50		0,00	0,00	0,00	2,26	0,00	0,00	2,26	10
27.12.2011 08:41:51		0,00	0,00	0,00	4,87	0,00	0,00	4,87	10
27.12.2011 08:41:57		0,11	0,00	2,64	0,46	0,00	0,00	3,21	10
27.12.2011 08:41:58		0,24	0,00	2,67	0,13	0,00	0,27	3,04	10
27.12.2011 08:41:59		1,34	0,00	2,64	0,00	0,00	0,00	3,98	10
27.12.2011 08:41:59		0,00	0,68	3,78	0,00	0,00	0,09	4,46	10
27.12.2011 08:41:59		1,23	0,00	1,60	1,23	0,00	0,00	4,05	10
27.12.2011 08:42:00		2,90	0,00	1,36	0,00	0,00	0,09	4,26	10
27.12.2011 08:42:00		0,95	0,00	2,92	0,14	0,00	0,14	4,01	10
27.12.2011 08:42:01		0,38	0,00	3,08	0,00	0,00	0,51	3,46	10
27.12.2011 08:42:01		0,48	0,00	3,57	0,00	0,00	0,00	4,05	10
27.12.2011 08:42:01		1,63	0,00	2,79	0,00	0,00	0,00	4,43	10
27.12.2011 08:42:02		0,71	0,00	3,36	0,36	0,00	0,00	4,43	10
27.12.2011 08:42:02		1,51	0,00	3,02	0,00	0,00	0,00	4,54	10
27.12.2011 08:42:03		1,53	0,00	2,92	0,00	0,00	0,21	4,45	10


Delete Data – advanced functionality only

To delete data, mark the datasets as described above.


Go to “EDIT → DELETE DATA” or use right click on the marked datasets:

0,07	1,59
0,00	1,66
0,00	2,37
History of Parameters	1,69
Edit Comment	0,00
Delete Data	0,25
Data Fields Selection	0,62
0,00	3,86
0,00	0,00
0,00	0,00
0,11	0,00

Click OK to delete the data.



## PARAMETERS MENU



Item	Function	Available for
History	Display parameters of datasets stored in the database	All instruments
Template	Apply a user-defined selection of parameters to an instrument	All but AlgaeTorch and BenthosTorch
Online	Display and change parameters in the instruments	All but AlgaeTorch and BenthosTorch
Recalc	Recalculate the datasets in the database with another set of calibration parameters	All but AlgaeTorch and BenthosTorch
Parameterset	Select a parameter set from the ones stored in the instrument – see description below.	FluoroProbe with parameter sets

The first 4 selections in the menu lead to one tab of the Parameters window.


**Please note:**

**Only parameters of the current type of instrument are shown.**


**It depends on the user level, whether a parameter is shown or not. The higher the user level, the more parameters are shown.**



For example of parameter listings for user levels, please see the examples below.

Access level “advanced functionality”:



Access level “User”:



Depending on the parameters and on the access level, the parameters are read only  or read and write .

## Overview

Recalculation with changed parameters.  
Create an editable parameter set by using "Create copy".  
Start recalculation by "Recalc".  
Restore original results by using "Undo".  
Calculation will effect all data currently in the view.

Creation date of parameter set 02.01.2012 11:20:27 (Test) **Adapted**

Name	Value	Unit
✓ Warm-up time	5	s
✓ Measurement interval	0	s
✓ Number of measurements	1	
✓ LED measuring time	20	1/10s
✓ Measurement time	20	1/10s
✓ Send results	on	
✓ Power-down after measurement	off	
✓ Store results	on	
✓ Start measurement on power-up	off	
✓ Calculate concentrations after each measurement	off	
⊖ Serial number	1709	
⊖ Software version	227	

### Type of operation

The selection of the type of operation corresponds to the commands in the menu. Depending on the selection, different types of operations can be done.

### Explanations

Gives a short introduction to what can be done.

### Buttons

Export: exports the current parameter set to another database for recalculation

Clipboard: copies the current parameter set to the clipboard to paste it into another application such as Office Word or Excel etc.

Upload: sends the parameters to the connected instrument

Matrix: shows all the calibration parameters of the connected instrument in one window.

### Name of parameter set

Each set of parameters can be identified by its time and date. Edited parameter sets or templates have an additional name given by the user.



### Type of parameter

Each tab shows the parameters of a specific topic.

### Values

Shows the current values of the parameters. Depending on the type of operation the value can be edited.

### Type of parameter set

Indicates the type of parameter set:

Original: this parameter set has been used originally for a measurement, it cannot be deleted.


Adapted: this is a parameter set has been adapted by the user, it can be used to recalculate the datasets.

Shows the current values of the parameters. Depending on the type of operation the value can be edited.

### History – advanced functionality only


Displays the parameters of the active dataset.

The history is also available on the Context Menu in table view:




The screenshot shows a table window titled "New.bdb [FP-17-09] :1". The table has three columns: "Date/Time [date]", "Comment", and "Green". The data rows contain timestamps from 27.12.2011 08:41:23 to 27.12.2011 08:42:00. A context menu is open over the row with timestamp "27.12.2011 08:41:24". The menu items are: "History of Parameters", "Edit Comment", "Delete Data", and "Data Fields Selection".

Date/Time [date]	Comment	Green
27.12.2011 08:41:23		
27.12.2011 08:41:24		
27.12.2011 08:41:24		
27.12.2011 08:41:24		
27.12.2011 08:41:24		
27.12.2011 08:41:24		
27.12.2011 08:41:59		
27.12.2011 08:41:59		
27.12.2011 08:42:00		
27.12.2011 08:42:00		



All parameter sets in the database can be displayed. To choose one parameter set, select the date from the drop-down box:




### Template Button

The „Template“ button is used to generate a new template from the current parameter set using different schemes. A scheme is a predefined selection of parameter – for example all the fingerprints. The parameter of the selection will be copied to a new template. Parameters that are not in the scheme will not be copied.


This is more convenient, than entering the parameter to a new template by hand.

The new template can be edited as any template.


First step is to choose the scheme with the selection of parameters.



Giving the new template a name:




Afterwards the new template can be displayed and edited. It can be used as any other template.




## Templates

The purpose of parameter templates is to store a partial set of parameters (for example offsets or measuring times or air pressure handling) in one set. This makes it easy to apply different setting to an instrument, depending on the location.




### Create a new template

Click new, all parameters of the current user-level are shown and empty. To change an entry, double click on the value box. When entering the first value, a name for the parameter set is requested:



The name and date is shown in the headline of the table:



Depending on the type of parameter, either an input box or a drop-down box is used:

<input checked="" type="checkbox"/>	Name of class b		
<input checked="" type="checkbox"/>	Name of class 7		
<input checked="" type="checkbox"/>	LED measurement interval	<input type="text"/>	\$
<input checked="" type="checkbox"/>	Enabled for fit flag class 1		
<input checked="" type="checkbox"/>	Enabled for fit flag class 5		
<input checked="" type="checkbox"/>	Enabled for fit flag class 6	<input type="text"/>	
<input checked="" type="checkbox"/>	Enabled for fit flag class 7	<input type="text"/>	off
<input checked="" type="checkbox"/>	Are Yellow Substances class 0		on
<input checked="" type="checkbox"/>	Are Yellow Substances class 1		


Click OK to store the template.

### Upload a template


1. connect the instrument to the PC ("TOOLS → CONNECT COM PORT")
2. go to "PARAMETER → TEMPLATES"
3. select the template from the drop-down box:

4. Click "UPLOAD":

bbe++ loads all the parameters from the instrument:



and sends the changed parameters back to the instrument.



All parameters not defined in the template will not be overwritten in the instrument.

### Delete a template

- select the template to be deleted
- click delete


## Online

The section “Online” is used to change the parameters of the instrument.


This item is available as soon as an instrument is connected.

1. connect the instrument to the PC (“TOOLS → CONNECT COM PORT”)
2. go to “PARAMETERS → ONLINE”

bbe++ loads all the parameters from the instrument:



The parameters are shown:




3. change the parameters as described in the “Template” section.
4. click OK to upload the new parameters
5. Use the “clipboard” button to copy all the parameters to the clipboard and paste it into EXCEL for example:

	A	B	C	D	E	F
1	Date: 09.11.2011 15:52:32					
2	Device: FP-17-09					
3						
4						
5	Process parameters					
6						
7	Name	Value	Unit			
8	Enter comment at start of measurement	on				
9	Air pressure from first measurement	off				
10						
11						
12	Common parameters					
13						
14	Name	Value	Unit			
15	Measurement interval	0	s			
16	LED measuring time	20	1/10s			
17	Measurement time	20	1/10s			
18	Power-down after measurement	off				
19	Start measurement on power-up	off				
20	Serial number	1709				
21	Software version	227				
22						
23						
24	Fit parameters					
25						
26	Name	Value	Unit			
27	Enabled for fit flag Green Algae	on				
28	Enabled for fit flag Bluegreen	on				
29	Enabled for fit flag Diatoms	on				
30	Enabled for fit flag Cryptophyta	on				
31	Enabled for fit flag sdsds#4	off				
32	Enabled for fit flag #5	off				
33	Enabled for fit flag #6	off				
34	Enabled for fit flag Yellow substances	on				
35	Date of calibration	21.11.2007				
36						
37						
38	Turbidity parameters					
39						
40	Name	Value	Unit			
41	Turbidity compensation	off				



## Recalc – advanced functionality only


The recalculation option is available in all instruments except for the AlgaeTorch. For the BenthosTorch the recalculation is limited to the recalculation by using another parameter set from the BenthosTorch. Editing of the parameter sets and templates are not available for the BenthosTorch.





To change the settings and recalculate the results, the following steps have to be executed:

1. Select the instrument and time period ("WINDOW → SELECT DATA").
2. Go to "PARAMETERS → RECALC"
3. Select a new dataset from the dropdown box or create a new one by editing a copy. Please note: only copies of an original parameter set can be edited. These sets are marked with "adapted".






4. To create a new parameter set (not in the BenthosTorch), select one, click “CREATE COPY” and enter a name for this parameter set.



5. Select the new parameter set and edit the settings as described in the “Template” section. Please note the number and type of parameters shown depends on the user-level. Only the parameters marked with the green symbol can be changed.



6. Click “APPLY” to store the changes in the new parameter set.
7. Click “RECALC” to recalculate the data with the new parameters.
8. Recalculated data sets are marked in table view at the beginning of each line:

Date/Time [date]	Comment	Gr
12.10.2011 11:25:14		0
12.10.2011 11:25:36		
12.10.2011 11:25:58		

### Restore data with original parameters


The original data of the measurement is not deleted by the recalculation. It is always possible to see the results calculated with the parameters valid when the measurement was taken.

To restore the original data:

1. Select the instrument and time period (“WINDOW → SELECT DATA”)
2. Go to “PARAMETERS → RECALC”
3. Click “UNDO”. Please note: This item is only available when recalculated data are available in the selected time period.




## VIEW MENU



### Toolbar

Displays or hides the toolbar.

The toolbar allows quick access to some important instructions:



### Status Bar

Displays or hides the status bar.




### Measurement

Displays / hides the measurement toolbar.


The toolbar "Measurement" is different for each type of instrument. It allows quick access to functions such as starting and stopping the measurement or retrieving the parameters.

The measurement toolbar of the FluoroProbe is given here as an example:



## Macro

The toolbar “Macro” contains buttons to access to the user-defined view macros. For details, please see the chapter “TOOLS” → “SETTINGS” → “MACRO MANAGER”.



For the grey buttons, no macro has been defined.

## Data Navigation

The toolbar “Data navigation” contains buttons to navigate within the data. “One page” corresponds to the time scale of the Concentration vs. Time graphs.



- ◆ go to the first dataset
- ◆ one page backwards
- ◆ one page forwards
- ◆ go to the last dataset

## Period

The toolbar “Period” contains buttons to scale the time axes to a given value:



- ◆ 1 hour
- ◆ 6 hours
- ◆ 1 day
- ◆ 7 days

## Session

The toolbar “Session” contains buttons to navigate from one measurement to another. One measurement is defined as a number of consecutive data sets that do not differ in the point of measurement by more than a given value. This value can be set in “TOOLS” → “SETTINGS” → “DISPLAY”.



## Connection State

The toolbar “Connection state” shows the status of the connection to the connected instrument:



The colour indicates the status:

- |                   |                                                               |
|-------------------|---------------------------------------------------------------|
| Black:            | Not connected, the PC has not tried to establish a connection |
| Yellow, blinking: | Not connected, the PC is trying to establish a connection     |
| Green:            | Connected                                                     |

To open the COM port setting click the connection state button

### Auto-Update

If auto-update is activated and a new data-set is generated during the measurement, this data-set is shown in the graphic or list window. If the displayed period does not contain this current point in time, the displayed period is changed automatically, so that the latest data-set is shown.

If auto-update is deactivated, the displayed period does not change. This is useful when doing data analysis during the measurement.

### Application look

Changes the style of the application


### Window – Advanced functionality only

This menu contains commands to show, arrange and edit different views to display the data.



### New Table Window

Opens a new window and displays the data according to the selected window type.




(See below)

Please see the list window “Concentration” of an AlgaeTorch as an example.

Date/Time [date]	Cyano [µg/l]	Conc Total [µg/l]	Turbidity [FTU]	Depth [
07.04.2010 11:33:33	36,1	77,5	0	
07.04.2010 11:33:47	35,8	76,9	0	
07.04.2010 11:34:02	35,8	76,9	0	


### New XYGraph Window

Opens a new window and displays the data according to the selected window type.




There are some view types pre-defined for each instrument. It is also possible to add new window types by using the view editor.

The following diagram shows the development of different algae classes measured with the bbe FluoroProbe.



### New Single Value View

Opens a new window and displays the data of one data-set according to the selected window type.



The data-set is shown in larger digits to help the user read the display even in bright sunlight. The latest data-set is shown.

Date/Time [date]	Green Algae [ $\mu\text{g/l}$ ]	Bluegreen [ $\mu\text{g/l}$ ]
22.08.2012 10:41:51	0,00	0,00
Diatoms [ $\mu\text{g/l}$ ]	Cryptophyta [ $\mu\text{g/l}$ ]	Yellow substances [ $\mu\text{g/l}$ ]
0,32	1,59	1,32
Total conc. [ $\mu\text{g/l}$ ]	Transmission [%]	Depth [m]
1,91	100,00	38,27
Temp. Sample [ $^{\circ}\text{C}$ ]		
6,15		

The “Single Value View” can be adapted to the screen by changing size and/or proportion of the window:

Date/Time [date]	Green Algae [ $\mu\text{g/l}$ ]
22.08.2012 10:41:51	0,00
Bluegreen [ $\mu\text{g/l}$ ]	Diatoms [ $\mu\text{g/l}$ ]
0,00	0,32
Cryptophyta [ $\mu\text{g/l}$ ]	Yellow substances [ $\mu\text{g/l}$ ]
1,59	1,32
Total conc. [ $\mu\text{g/l}$ ]	Transmission [%]
1,91	100,00
Depth [m]	Temp. Sample [ $^{\circ}\text{C}$ ]
38,27	6,15

Date/Time [date] 22.08.2012 10:41:51
Green Algae [ $\mu\text{g/l}$ ] 0,00
Bluegreen [ $\mu\text{g/l}$ ] 0,00
Diatoms [ $\mu\text{g/l}$ ] 0,32
Cryptophyta [ $\mu\text{g/l}$ ] 1,59
Yellow substances [ $\mu\text{g/l}$ ] 1,32
Total conc. [ $\mu\text{g/l}$ ] 1,91
Transmission [%] 100,00
Depth [m] 38,27
Temp. Sample [ $^{\circ}\text{C}$ ] 6,15

Date/Time [date]	Green Algae [ $\mu\text{g/l}$ ]	Bluegreen [ $\mu\text{g/l}$ ]	Diatoms [ $\mu\text{g/l}$ ]	Cryptophyta [ $\mu\text{g/l}$ ]
22.08.2012 10:41:51	0,00	0,00	0,32	1,59
Yellow substances [ $\mu\text{g/l}$ ]	Total conc. [ $\mu\text{g/l}$ ]	Transmission [%]	Depth [m]	Temp. Sample [ $^{\circ}\text{C}$ ]
1,32	1,91	100,00	38,27	6,15

### Cascade

Arranges the windows one behind the other.




## Tile Horizontally/Vertically

Initiates the organization of the screen into mutually non-overlapping frames.

## Select Data

The “Select Data” window can be used to select data of a specific instrument and a specific time range from the current database.




Type: shows the types of instruments available in this database.

Unit: shows the serial numbers of the type of instruments chosen above available in this database.

## Time range

Selects the time range of data to be selected. It is possible to select a fixed time range such as “last month” or to do a manual selection by entering a range FROM “date” TO “date”.

***Hint: Further selections can also easily be made in a graphic window.***



After selecting the time range, the data is shown in the default view.

Date/Time [date]	Conc 0 [Green Algae] [µg/l]	Conc 1 [Bluegreen] [µg/l]	Conc 2 [Diatoms] [µg/l]
22.02.2007 15:27:25	0	1,34	
22.02.2007 15:27:34	0	1,43	
22.02.2007 15:27:40	0	1,41	
22.02.2007 15:27:47	0	1,37	
22.02.2007 15:27:53	0	1,45	
22.02.2007 15:28:00	0	1,4	
22.02.2007 15:28:06	0	1,34	
22.02.2007 15:28:13	0	1,35	
22.02.2007 15:28:20	0	1,38	
22.02.2007 15:28:26	0	1,42	
22.02.2007 15:28:35	0	1,3	
22.02.2007 15:28:41	0	1,47	
22.02.2007 15:28:48	0	1,4	
22.02.2007 15:28:54	0	1,33	
22.02.2007 15:29:01	0	1,46	
22.02.2007 15:29:08	0	1,39	
22.02.2007 15:29:14	0	1,37	

### View Editor

The View Editor is used to generate your own table or graphic window type as well as a template for the export to text files. The View Editor is an assistant that leads through the whole process.

When starting the View Editor, the current view is selected automatically to be edited or used as a template for a new view. Other views may be selected.


For each instrument there are pre-defined views. These views can be edited but not deleted. Changes in pre-defined views can be reset to their original status. User-defined views can be edited and deleted.

### Generating and editing a table view window

To explain the function of the View Editor, in the following steps a table view showing date/time, total chlorophyll concentration, Bluegreen concentration and depth is generated as an example.

### View Selection

1. Open the view editor:




2. Choose the type of instrument the view is to be used for (in the example "FluoroProbe").


3. Choose the type of view that is to be generated or edited (in the example “Table view / single value view”).
4. Choose one existing view to edit or as a template for a new one.  
In case of a pre-defined view, the view can be reset to the original status (“RESET”).  
In case of a user-defined view, the view can be deleted (choose a custom view and click “DELETE”).
5. Click “NEXT”.

### Field Configuration

This window shows all available data of the selected instrument. The left column shows the data that is currently not used in the chosen view. In the right column are data that will be displayed in the resulting view. The order in the right column corresponds to the order in the table.



6. Click on “Green Algae [µg/l]” and “← REMOVE” to remove the concentration of green algae from the view.
7. Repeat with all entries but “Date/Time (date)”; “Bluegreen [µg/l]”, “Total conc. [µg/l]” and “Depth [m]”.
8. Click on “Total conc. [µg/l]” and ”UP” to change the order.




9. Click "Next".

### General Properties


The General Properties page allows the user to change the background color of a "Single value" view.




10. Click on the color picker box behind „Background color“ and choose a new color.



The chosen color is now shown in the color picker box




Example of a change in background color:



### Save View


11. Enter a name for the view.



12. Click “OK” to confirm. The current data is then shown with the new view.

Date/Time [date]	Total conc. [µg/l]	Bluegreen [µg/l]	Depth [m]
22.08.2012 10:40:18	3,21	0,00	0,82
22.08.2012 10:40:18	3,04	0,00	1,12
22.08.2012 10:40:20	3,98	0,00	1,74
22.08.2012 10:40:20	4,46	0,68	2,14
22.08.2012 10:40:21	4,05	0,00	2,45
22.08.2012 10:40:22	4,26	0,00	2,76
22.08.2012 10:40:23	4,01	0,00	3,16
22.08.2012 10:40:23	3,46	0,00	3,47
22.08.2012 10:40:24	4,05	0,00	3,78
22.08.2012 10:40:25	4,43	0,00	4,08
22.08.2012 10:40:26	4,43	0,00	4,39
22.08.2012 10:40:26	4,54	0,00	4,69
22.08.2012 10:40:27	4,45	0,00	5,00
22.08.2012 10:40:28	4,30	0,00	5,31
22.08.2012 10:40:29	4,33	0,00	5,61
22.08.2012 10:40:29	4,22	0,00	5,92
22.08.2012 10:40:30	4,34	0,00	6,23
22.08.2012 10:40:31	4,03	0,00	6,53
22.08.2012 10:40:32	4,29	0,00	6,74
22.08.2012 10:40:33	4,14	0,00	7,14

13. The new view is now available in the Window menu.





### Generating and editing a graphic view window

To explain this function, in the following steps a graphic view showing date/time and total chlorophyll concentration is generated as an example.

#### View Selection


1. Open the view editor:




2. Choose the type of instrument the view is to be used for (in the example “FluoroProbe”).
3. Choose the type of view that is to be generated or edited (in the example “Graphics view”).
4. Choose one existing view to edit or as a template for a new one (in the example “Conc. / depth”).  
In case of a pre-defined view, the view can be reset to the original status (“RESET”).  
In case of a user-defined view, the custom view can be deleted (“DELETE”).
5. Click “NEXT”.

Independent Axis

6. Choose the design according to the given examples (in the example Y-Axis)


Field Configuration

7. Select the data for the axis chosen and click on “ADD→”. Only one of the available entries is possible (in the example “Depth”)






8. After choosing one entry all other entries vanish. Removing the entry by clicking “←REMOVE” brings back all available entries.



9. Click on „PROPERTIES“

#### Axis Properties - first axis

Some properties can be set for each axis. For the first axis, this is the position and the direction.





Normally the first axis is on the lower or left side of the graph. The direction is normally defined as:


X-Axis: left side → low values                      right side → high values

Y-Axis: bottom → low values                          top → high values


To display e.g. chlorophyll concentration versus depth, it is easier to use the top position of the X-axis and to invert the direction of the Y-axis for the depth. Please see the following example: the first picture shows the standard settings, the second one with a direction-inverted Y-axis and inverted position of the X-axis.




10. Click "OK" to confirm the settings.
11. Click "NEXT".
12. Now the data and settings for the other axis can be selected.

*Field Configuration - second axis*

13. After selecting one type of data, only those types which have the same unit(s) remain in the right column.




14. For example, "Total chlorophyll concentration" and "Bluegreen" have to be selected here.



### Axis Properties - second axis

For the second axis, the properties can be set as shown above. Additionally available scaling options are:



**Manual scaling:** enter fixed values for the minimum and the maximum of this axis - enter 0 for both to disable this feature.

**Auto-scale on:** enables the auto-scaling feature - the axis will be scaled so that all selected data are visible.

**Auto-scale limits:** these settings are used to obtain a nice looking graph in two special situations:

- datasets with very low noise:

In the auto-scaling mode, the graph looks as if there are extreme variations, but the range is very small due to the auto-scaling. In this case, it is better to use the option:

lower limit: lower than and

upper limit: higher than

This leads to a minimum span in the graph.

- datasets with outliers:

In this case, it is difficult to analyze the data because the auto-scaling generates a high span value. To prevent this, please use:

lower limit: higher than for outliers with low values and

upper limit: lower than for outliers with high values


This leads to an optimized span in the graph.

15. Click "Next"

### Attributes

This sets the color and size of the dots and lines as well as the color of the background.

16. Set the value in the attributes window according to your needs.




Pixel number: enter the dot size of the value points.


Background color: click on the colored button to change the background color of the graph.

For illustration, please see the examples with connected and unconnected dots, different sizes and background colors.

Different dot sizes:



Different background colors:




### Column Attributes

Select a dataset.


Color: select a dataset and click on the colored button to change the color of the dots.

Connected: check the box to connect the dots in the graph.




**Hint: The data points are not connected if the time of measurement differs more than the given value in "TOOLS → SETTINGS → DISPLAY".**

Different dot colors:




Unconnected dots / connected dots:




17. Click "Next"

Save view

18. Enter a new name for the view. It is not possible to edit pre-defined views. Please store a view based on a pre-defined view with a new name.




19. The current data is shown using the new view after clicking "OK".



## Generating and editing an ASCII Export

### View Selection


1. Open the view editor:



2. Choose the type of instrument the export is to be used for (in the example “AlgaeTorch”)
3. Choose the type of export that is to be generated or edited (in the example “ASCII export”)
4. Choose one existing view to edit or as a template for a new view. The delete button is used to delete a view. It is not possible to delete pre-defined views.
5. Click “NEXT”



### Field Configuration

This window shows all available data of the selected instrument. The left column shows data which are not currently used in the chosen view. In the right column are data that are to be displayed in the view. The order in the right column corresponds to the order in the table.



6. Select the required datasets from the right column. Use the up and down buttons to determine the order of the data in the export table.
7. Click “Next”.






8. Define separator and end of line characters according to the needs of the importing software.
9. Click "Next".

### Save View

10. Enter a new name for the view. It is not possible to edit pre-defined views. Please store a view based on a pre-defined view with a new name.




11. The new view is now available in the "File" menu.



### Comment Fields

For each type of instrument, special comments can be defined. These comment fields can be used at the beginning of the measurement. The comments names and types can be defined within this dialog.



Type of instrument: the type of instrument is defined in the settings "TOOLS → SETTINGS → ACCESS".


Used:	if tagged, the comment can be used at the start of the measurement, otherwise the comment is hidden.
Description:	name of the comment
Type:	“Text” allows entry of text, “Number” allows digits only; this is relevant for the way the comments are sorted in a list view.
Template:	stores the comment definition to transfer it to a new database; the transfer is done when opening a new database.

### 1. ... 2. ... 3. ...


Shows all the currently opened windows.

## TOOLS

Menu (advanced functionality)



Menu (basic functionality)




## Settings

Here, all the preferences for the general functionality of bbe++ can be set.

## Database

The 'Database' tab shows the currently used database. The checkbox can be used to open this database after starting bbe++. We recommend you to store all the data in the same database.



### Automatic start of measurement

For continuously working instruments (AOA or FluoroProbe / AlgaeTorch) the software can be configured to start the measurement directly after starting the software. The feature can be used to ensure, that the measurement starts again after a power fail.


### Auto Archive

If there are many datasets in the database, navigation in the database becomes slow. The "Auto Archive" function automatically archives older data in the database. If there are data older than – in the example – 60 days, the auto-archive function is started. Only datasets that are a maximum of 30 days old remain in the database.

The archived data is written in a database similar to an exported database. The name corresponds to the date of the first dataset in this database.

## Display


Enter the time period which defines the maximum interval between two measurements of the same series. The setting is used for the “Next measurement” button in the toolbar.




Next and previous measurement buttons in the toolbar.



Furthermore, this setting is used to control the drawing of lines between two dots. If the interval between two measurements is shorter than the given interval, a line is drawn – otherwise not:






Minimal interval set to 1 minute

### Access

Depending on the access level selected, the software shows more or fewer options. The main difference is the number of items shown in the menus and the number of parameters shown.

If an instrument with editable parameters is connected, the same setting is maintained for editing.




User-levels:

Level	Access
Basic functionality	<ul style="list-style-type: none"><li>• Just basic operation of the instrument</li><li>• Measurements with pre-defined parameters only</li></ul>
Advanced functionality	<ul style="list-style-type: none"><li>• All menus accessible / all parameters customizable</li></ul>

The “Advanced functionality” level can be password-protected. As soon as the checkbox is ticked, a password can be entered. Afterwards the password is required to enter the Advanced functionality level.

***Please note:***

***The “Advanced functionality” level remains until “Basic functionality” level is chosen again. Even when restarting the software, the “Advanced functionality” level and a given password will still be valid.***

Default Device


The bbe++ software is able to control different types bbe instruments, but only one type at one time. The type of instrument is entered here. Various settings within the bbe++ software can only be made for the default type instrument. To make entries for an other type of instrument, the default type has to be changed before.

Offline Mode:

Choose the offline mode to display data of different type of instruments. If the offline mode is active, no measurement is possible.


## Macro manager

The bbe++ software provides the option to store the arrangement of the view windows. One of the stored views can be used as the default view automatically displayed after starting bbe++.



To generate a new view macro:


1. Arrange the windows.
2. Enter the name of the new view.
3. Check the "Use at program start" box as required.
4. Confirm.



**Logger settings**

The logger settings are used in case of problems with the bbe++ software. Please ask the bbe service for advice.


Depending on the checked options more or less information is stored in the file “error.log” in the program folder of bbe++ software.



**Shortcuts**

With the F8 key, one single view can be retrieved. Please enter the single view macro to be used for this.






Date/Time [date] 26.06.2007 12:19:31	Green Algae [ $\mu\text{g/l}$ ] 0,00
Bluegreen [ $\mu\text{g/l}$ ] 0,77	Diatoms [ $\mu\text{g/l}$ ] 0,00
Cryptophyta [ $\mu\text{g/l}$ ] 0,00	Yellow substances [ $\mu\text{g/l}$ ] 0,00
Total conc. [ $\mu\text{g/l}$ ] 0,77	Transmission [%] 98,22
Depth [m] 0,11	Temp. Sample [ $^{\circ}\text{C}$ ] 21,61

F8 shows the selected macro


### Language

This tab shows a list of the available languages of the bbe++ software. To change the language, choose the language you wish and click "OK". The bbe++ software must be restarted for the changes to take effect.



### IO moduls

This item is used to connect an IO module (for example 4-20mA or relay outputs) to the PC.




- the COM port can either be chosen from the drop down box or the software can search all COM ports for the selected instrument.
- once an instrument has been found the COM port is stored for the next use.
- To configure the devices goto Tools -> device -> ... . This item is not available for all instruments.

With "search MCU" a multi channel unit for the AlgaeOnlineAnalyser can be connected.

## Text results


This item is used to do the configuration of writing text files and sending data via RS232.



## File output

If the file output is configured daily text files are written by bbe++. The files are stored in the path given in “file path”.

The name is the serial number of the instrument date:



The format of the output is defined by a “view” of the type “ASCII export”. For details please see the chapter “view editor”. In the view editor the data columns and the format can be defined.

Files older than given in the “max. age of files” are deleted automatically. If this is set to “0” no files will be deleted.

## Enable advanced/basic functionality

Change between advanced and basic functionality. See also Tools -> Settings -> Access.

## Connect COM port



This item is used to connect an instrument to the PC.

**Please note:**

- bbe++ is looking for an instrument of the selected type only (Tools → Settings → Access).
- the COM port can either be chosen from the drop down box or the software can search all COM ports for the selected instrument.
- once an instrument has been found the COM port is stored for the next use.

### Device

The “Device” submenu varies depending on the selected type of instrument. Each type has its own submenu.


### DemoMode

To switch the bbe++ software to the demo mode select FluoroProbe (Tools -> Settings -> Access) and disconnect all instruments. Click on “DEMO mode” in the device submenu. Now, the software can be operated as if a FluoroProbe were connected. Some procedures such as calibrations are limited in function.


To reset the demo mode, please restart bbe++.



Once the DEMO mode is selected, the submenu changes to the FluoroProbe entries:



### AOA device Submenu



#### Enter comment at start of measurement:

Check to enter one or more comments when starting the measurement procedure.

#### Activity:


Check to perform activity measurements. If unchecked the activity measurement will be skipped.

Continuous:

Check to perform continuous measurements (standard). If unchecked, only one measurement will be done. This can be used for lab operation.

Test cleaning:

A cleaning procedure can be triggered.



Up/down:

Starts the cleaning piston to move up or down. The travelling time is indicated (normally around 20s)

Test pump:

The pump can be operated manually.



The test pump dialogue can also be used to check the required pumping time to get fresh sample to the instrument.

Analog output (option):


The analog outputs can be configured according to individual needs. The outputs can be used for different parameters and the range for the results can be set for 4-20 mA.

Additionally the sample source can be filtered:

None: All data of the chosen parameter will be applied, in depended from the sample source

Source n: Only data of the sample source n will be applied.

Depending on the number of installed analogue outputs more or less tabs will be displayed. In the screenshot below 8 analogue outputs are installed.



In this example, on channel 1, 4 mA is supplied at a chlorophyll concentration of 0 µg/l and 20mA at a concentration of 100 µg/l – for any sample source. Values that are above or underneath these are converted to 20 respectively 4mA.

On channel 2, 4 mA is supplied at a chlorophyll concentration of 0 µg/l and 20mA at a concentration of 200 µg/l of green algae – for sample source 2.

The current output of channel 1 will only be refreshed after the measurement of sample source 1. The current output of channel 2 will only be refreshed after the measurement of sample source 2.


Example for the setting above:

Source measured	Conc. Source 1 [µg/l]	Conc. Source 2 [µg/l]	Output channel 1 [µg/l]	Output channel 2 [µg/l]
			No data filter	Source 2 only
2	50	0	4mA	4mA
1	50	0	12mA	4mA
2	50	100	20mA	10mA
1	50	100	12mA	10mA
2	0	100	20mA	10mA
1	0	100	4mA	10mA

There are also output options for the alarms. If there is no alarm, the output is set to 4mA, otherwise to 20mA.

The Test 4 mA, Test 12 mA and Test 20 mA buttons can be used for test purposes. If one of these buttons is clicked, the current indicated on the button can be measured at the corresponding output.


Depending on the output device, the 4 to 20mA output may also be calibrated. If this option is available, just click on the 4 and 20mA buttons in the calibration section of each output:



Using the arrow buttons it is possible to adjust the corresponding mA value so that the SCADA system shows the same value. Press OK to confirm the settings.

### Relay Output (Option)

The relay outputs can be configured according to individual needs. The menu item is only available if a relay board is installed.



The relay outputs can be set in the following modes:

Disabled	No function
WatchDog	The relay changes between on and off one time per minute during the measurement. No toggling if a hardware alarm occurs.
Conc. alarm of algae class...	The relay switches to the ON position if an alarm of the specified algae class occurs.
Conc. alarm of total chlorophyll	The relay switches to the ON position if an alarm of the total chlorophyll concentration occurs.
Any conc. alarm	The relay switches to the ON position if any concentration alarm occurs.
Hardware alarm	The relay switches to the ON position if any hardware alarm occurs.







	525	570	610	590	470	UV LED
Old offset	0,068964	0,21124E	0,15214E	0,21737E	0,114701	0,45456E
New offset	0	0	0	0	0	0

- Click "START"

Please insert a cuvette with ultrafiltrated sample water into the device, then press OK.

- Fill a clean cuvette with distilled water, place it in the AlgaeLabAnalyser and confirm with "OK".
- Now 10 measurements of the sample will be executed. The number of the current measurement will be shown in the field "Measurement"
- The new gradient for each LED will be shown. "APPLY" starts the transfer of the new parameters to the AlgaeLabAnalyser.

	525	570	610	590	470	UV LED
Old offset	0,08067	0,22131	0,19676	0,27305	0,13140	0,49071
New offset	0	0	0	0	0	0

"MORE" shows the settings for the digital-analog converter that controls the LED. This information is only of value for a maintenance technician.

### Chlorophyll-A

This menu item enables an offset measurement with the aid of a filtered sample. In comparison to the other offset measurements, this is characterized by a significantly faster uptake of the offset, since the F0 and Fm measurements are excluded. Please use this offset calibration when the yellow substance measurement is excluded.

### Offsets

The offset is the signal without algae that has to be subtracted during the measurement. Two offsets have to be calibrated:

1. The offset in distilled water. This water has to be free of algae **and** free of yellow substances. This offset is used when the determination of yellow substances is switched on. The concentration of yellow substances is determined during the measurement.
2. The offset in ultrafiltrated water. This water has to be free of algae. This offset is used when the determination of yellow substances is switched off. The concentration of yellow substances is “in” the offset and so filtrated sample water has to be used for the offset.

### Calibration solution for the „ultra filtrated water offset and the yellow substances

In the lab, the calibration of the yellow substances can be easily done with tap water. In cases of very small chlorophyll-a concentration, it is recommended to use ultrafiltered sample of the water (from the water to be examined). For this the following equipment is needed:

**Type of filter:** Nitrocellulose filter, pore size 0.45 µm, diameter 47mm (for example: Sartorius Stedim)

**Vessels:** Depending on the sample volume flask with 5 L content


Negative pressure of 800mbar. The volume of the sample is determined by the instrument to be calibrated, BBE. For an ALA be about 100ml, requires approximately 5L filtrate for an AOA / AlgaeGuard / A-Tox sensor / AlgaeTorch about 1L and for an immersion probe.

### Offsets (ultrafiltrated)

Performs a calibration of the ultrafiltrated sample water offset.

1. Go to “CALIBRATION → OFFSETS (ULTRAFILTRATED)”. The calibration parameters will be transferred from the AlgaeLabAnalyser and displayed.

2. Click on “START” to initiate the calibration



3. Please fill a clean cuvette with distilled water, place it in the AlgaeLabAnalyser and click “OK”.
4. Subsequently 3 x 10 measurements are executed:

- Offset of the chlorophyll concentration measurement (F-measurement)
- Offset of the so called fo and fm measurement (Genty parameter).

The current measurement is displayed in the “Type” field. The current number in the “Measurement” field.

5. The transmission is displayed for control purposes
6. After the calibration, the new calibration values are displayed. Please confirm with “Apply” to transfer the new calibration.

	525	570	610	590	470	UV LED	Trans.
Old offset	2,36589	1,13209	0,47290	0,45963	5,96178	5,21487	Digits
New offset	0,03124	0,05573	0,11227	-0,08894	0,00689	-0,15569	Digits %

“MORE” shows the settings for the digital-analog converter that controls the LED. This information is only of value for a maintenance technician.

### Offsets (distilled)

Performs a calibration of the distilled water offset.

1. Go to “CALIBRATION → OFFSETS (DISTILLED)”. The calibration parameters will be transferred from the AlgaeLabAnalyser and displayed.

	525	570	610	590	470	UV LED	Trans.
Old offset	0	0	0	0	0	0	digits
New offset	0	0	0	0	0	0	digits %

2. Click on “START” to initiate the calibration
3. Please fill a clean cuvette with distilled water, place it in the AlgaeLabAnalyser and click “OK”.

4. Subsequently 3 x 10 measurements are executed.:

- Offset of the chlorophyll concentration measurement (F-measurement)
- Offset of the so called fo and fm measurement (Genty parameter).

The current measurement is displayed in the “Type” field. The current number in the “Measurement” field.

5. The transmission is displayed for control purposes
6. After the calibration, the new calibration values are displayed. Please confirm with “Apply” to transfer the new calibration.

	525	570	610	590	470	UV LED	Trans.
Old offset	2,20743	1,49263	1,88671	0,84044	1,79128	0,22316	Digits
New offset	0,01745	0,08160	0,05635	0,02822	0,06357	-0,09353	Digits 95,37 %

“MORE” shows the settings for the digital-analog converter that controls the LED. This information is only of value for a maintenance technician.

## Fingerprints

This menu item subsumes the fingerprints of algae classes and yellow substances. The fingerprints contain the fluorescence characteristics of the algae classes or yellow substances. The fingerprint of the algae classes shows the signal of the pure solution based on a concentration of 1µg/l. For the fingerprints of yellow substances the concentration is chosen arbitrary. Their classification serves mainly the correction of the fluorescence signal during the classification of algae.

The names of the algae and if a fingerprint is an algae or a yellow substance can be defined in the parameters.


For the calibration of the algae fingerprints, please use pure solutions of the algae division with a known chlorophyll concentration. The concentration should be in the range of 20-80 µg/l.

1. Go to “CALIBRATION → FINGERPRINTS”. Now the parameters are imported from the device and the calibration dialogue is shown.

	525	570	610	590	470	UV LED	Trans.
Old factors	0,43353E	0,142361	0,19808E	0,13105E	1,03524E	1,06661E	
New factors	0	0	0	0	0	0	

2. Select the algae class to calibrate
3. Insert the concentration. Concentration values should be in a range of 20 – 80 µg/l.
4. If the concentration is determined later insert “1” and use the “calibration by factor” feature to adapt the calibration after specifying the concentration.


5. Please choose from the "Type of offset"-field if the sample was diluted with filtrated or distilled water.
6. Click on "START"



7. Please fill the calibration sample in a clean cuvette, place it in the AlgaeLabAnalyser and click "OK".
8. After the calibration the new fingerprint is shown.
9. Apply the values by clicking "Apply".
10. Repeat the procedure with all algae classes activated in the AlgaeLabAnalyser

### Correction factors

Go to "Calibration → Correction factors". Now the parameters are imported from the device and the fingerprint correction dialogue is shown.



### Correction of a fingerprint

This item allows the correction of a specific fingerprint using a particular factor. If one algae class is under- or overestimated compared to your reference, the fingerprint can be changed automatically. Entering for example the factor 2 will double the result for the chosen algae class. After entering the number, the new fingerprint is displayed. Use the "APPLY" button to send it to the sensor.

FACTORY DEFAULT: This option can be used to reset a specific fingerprint to the factory settings. This feature is only available in sensors with an internal software 1.96 or later. Otherwise, this option is greyed out.

### Calibration of Fingerprints with Subsequent Wet-Chemical Analysis:

If the AlgaeLabAnalyser is calibrated with a reference solution where the concentration is determined subsequently, the following steps have to be carried out:


- Do a fingerprint calibration as described above, but enter 1 for the concentration of the reference solution
- Do a wet-chemical analysis
- Use the “Calibration by factor” feature to adapt the calibration to the result of the subsequent wet-chemical analysis.

Example:

The result of the subsequent chlorophyll analysis is 35µg/l. In this case the factor is 35.


Attention: the AlgaeLabAnalyser is not ready for use until the last step has been carried out.

## HELP



## About ...

Displays the start-up screen.



For further questions, please contact our service. Please let us know the Software version you are using and the serial number of your instrument.

New version can be downloaded from <http://www.bbe-moldaenke.de/>.

## bbe++ manuals ...

List of currently installed bbe++ manuals. Depending on the selection made during the installation, there may be manuals in different languages and for different instruments.

To make further manual available, please store the PDF files in the following folder:

Windows XP:	C:\Programme\bbe++\
Windows Vista or Windows 7:	C:\Programme (x86)\bbe++\



## THE VIEWS

### TABLE VIEWS

A table view is opened by “WINDOW → NEW TABLE VIEW → TABLE VIEW NAME“.

Date/Time [date]	Conc 0 [Green Algae] [µg/l]	Conc 1 [Bluegreen] [µg/l]	Conc 2 [Diatoms] [µg/l]	Conc 3 [Cryptophyta] [µg/l]	Conc 4
15.06.2009 09:10:50	0	0,813	0	0	0,495
15.06.2009 09:10:53	0	1,04	0	0	0
15.06.2009 09:10:55	0	0	0	0	0
15.06.2009 09:10:57	0	2,28	0	0	0
15.06.2009 09:11:00	1,08	0,829	0	0	0
15.06.2009 09:11:02	0	0,142	0	0	2,09
15.06.2009 09:11:04	0,0475	0,956	0	0	0,908
15.06.2009 09:11:07	0	0,932	0	0	0
15.06.2009 09:11:09	0	1,2	0	0	0
15.06.2009 09:11:12	0	0,74	0	0	0
15.06.2009 09:11:14	0	0,561	0	0	0,936
15.06.2009 09:11:16	0	0,698	0	0	0
15.06.2009 09:11:19	0,945	0,779	0	0	0
15.06.2009 09:11:21	0	0,804	0	0	0
15.06.2009 09:11:23	3,35	0	0	0	0
15.06.2009 09:11:26	1,07	0,488	0	0	0
15.06.2009 09:11:28	0	0,261	0	0	2,03
15.06.2009 09:11:31	1,52	0,227	0	0	0
15.06.2009 09:11:33	0	1,1	0	0	0
15.06.2009 09:11:35	0	1,58	0	0	1,97
15.06.2009 09:11:38	5,41	0	0	0	0,883
15.06.2009 09:11:40	3,26	0	2,18	0	0
15.06.2009 09:11:42	0,792	0	2,82	0	0
15.06.2009 09:11:45	2,92	0	1,33	0	2,85
15.06.2009 09:11:47	6,68	0	0	0	1,54
15.06.2009 09:11:50	6,56	0	0,939	0	0,492
15.06.2009 09:11:52	0	0,111	1,16	0	0
15.06.2009 09:11:54	0	0,363	0	0	0,598

Within the list view, the right mouse button offers detailed information about the selected data set.

Date/Time [date]	Comment	Comment [1]
14.10.2015 14:53:58		
14.10.2015 15	History of Parameters	
14.10.2015 15	Edit Comment	
14.10.2015 15	Delete Data	
14.10.2015 16	Data Fields Selection	
14.10.2015 16		





## Delete data

To delete data rows from the data base, mark one or more row by clicking on the row.

Hold the Ctrl. key while clicking to add more rows to the selection.

14.10.2015 16:53:58
14.10.2015 17:03:58
14.10.2015 17:12:58
14.10.2015 17:53:58
14.10.2015 18:03:58
14.10.2015 18:12:58
14.10.2015 18:53:58
14.10.2015 19:03:58
14.10.2015 19:12:58
14.10.2015 19:53:58
14.10.2015 20:03:58
14.10.2015 20:12:58
14.10.2015 20:53:58
14.10.2015 21:03:58
14.10.2015 21:12:58
14.10.2015 21:53:58
14.10.2015 22:03:58

Hold the Shift. key while clicking to expand the.


14.10.2015 16:03:58
14.10.2015 16:12:58
14.10.2015 16:53:58
14.10.2015 17:03:58
14.10.2015 17:12:58
14.10.2015 17:53:58
14.10.2015 18:03:58
14.10.2015 18:12:58
14.10.2015 18:53:58
14.10.2015 19:03:58
14.10.2015 19:12:58
14.10.2015 19:53:58
14.10.2015 20:03:58
14.10.2015 20:12:58
14.10.2015 20:53:58
14.10.2015 21:03:58
14.10.2015 21:12:58
14.10.2015 21:53:58
14.10.2015 22:03:58
14.10.2015 22:12:58
14.10.2015 22:53:58
14.10.2015 23:03:58
14.10.2015 23:12:58
14.10.2015 23:53:58

Select "Delete Data"

14.10.2015 17:03:58
14.10.2015 17:12:58
14.10.2015 17:53:58
14.10.2015 18:03:58
14.10.2015 18:12:58
14.10.2015 18:53:58
14.10.2015 19:03:58
14.10.2015 19:12:58
14.10.2015 19:53:58
14.10.2015 20:03:58
14.10.2015 20:12:58
14.10.2015 20:53:58


History of Parameters
Edit Comment
Delete Data
Data Fields Selection

All data or the selected only can be deleted. Click OK to delete.



### Data Fields Selection

To add or remove columns from the list view use "Data fields selection". The dialog is similar to the dialog of the View Editor.



Please select or deselect the columns. The list view is updated immediately.

### Sorting of the list view

The list view is sorted by date/time. It can be sorted by any other column by clicking on the header of the column. Clicking twice leads to the reversed order.

Sorted by time/date ascending:

Date/Time [date]	Green Algae [µg/l]
26.06.2007 10:04:14	0,00
26.06.2007 10:04:29	0,00
26.06.2007 10:04:45	0,00
26.06.2007 10:05:00	0,00
26.06.2007 10:05:16	0,00
26.06.2007 10:05:32	9,48
26.06.2007 10:05:47	0,00
26.06.2007 10:06:03	0,00
26.06.2007 10:06:19	0,00

Sorted by time/date descending:

Date/Time [date]	Green Algae [ $\mu\text{g/l}$ ]
26.06.2007 12:19:30	0,00
26.06.2007 12:19:14	0,00
26.06.2007 12:18:59	0,00
26.06.2007 12:18:44	0,00
26.06.2007 12:18:29	0,00
26.06.2007 12:18:13	0,00
26.06.2007 12:17:58	0,00
26.06.2007 12:17:43	0,00
26.06.2007 12:17:28	0,00

Sorted by green algae descending:

Date/Time [date]	Green Algae [ $\mu\text{g/l}$ ]
26.06.2007 10:05:32	9,48
26.06.2007 10:09:42	2,97
26.06.2007 11:20:26	2,89
26.06.2007 11:28:16	1,98
26.06.2007 10:16:14	1,31
26.06.2007 11:08:26	1,24
26.06.2007 11:55:09	0,32
26.06.2007 10:17:48	0,13
26.06.2007 11:30:53	0,12


### XY GRAPH VIEWS

Features of the XY graph view:



1. Select data of interest
2. Scaling of the axes (right click)
3. Display of the corresponding parameters (right click)
4. Editing a XY view, to change colours the type of graph.
5. Navigation within the data

### Select data of interest


Please click and draw the mouse to select the required data.




A pop-up window is shown to choose the axis to which the selection is to be applied.




- X only                    only the selection on the horizontal axis is applied.
- Y only                    only the selection on the vertical axis is applied.
- X and Y                   the selection on the horizontal and the vertical axis is applied.



X only applied




Y only appliedMinimal



X and Y applied

### Right-click menu




#### Features:

- adjust the scaling .
- activate and deactivate of the autoscaling feature.
- show the parameters of a single data-set.
- return to the default scaling.

### Scaling of the axis

By default, scaling is done automatically. This means, all datasets of the selected period are displayed. To have a closer look at the data, it is sometimes useful to change the scaling. In the right-click menu, there is a “Scale axis” option:


In case of a time axis, the scaling tab of this axis is as follows:



The first and last point in time have to be entered.

**Hint: For more information on scaling an axis please see the "window" section in this manual.**

In case of a "non-time axis", more possible adjustments are available:



**Manual scaling:** enter fixed values for the minimum and the maximum of this axis - enter 0 for both to disable this feature.

**Auto scale limits** these settings are used to obtain a nice looking graph in two special situations:

- data sets with very low noise:

In the auto-scaling mode, the graph looks as if there are extreme variations, but the range is very small due to the auto-scaling. In this case, it is better to use the options:

lower limit:      lower than      and

upper limit:      higher than

This leads to a minimum span in the graph.

- data sets with outliers:

In this case, it is difficult to analyse the data because the auto scaling generates a high span value. To prevent this, please use the

lower limit:      higher than for outliers with low values      and

upper limit:      lower than for outliers with high values  
This leads to an optimised span in the graph.

### Autoscaling

Enables the auto-scaling feature. The axis will be scaled so that all selected data are visible.

### Parameters

Show the parameters of the data-set

### Default scaling

Click here to undo all the changes in the scaling options.




## PARAMETERS

The following pages show the parameters that are visible in the service access level. In the user level not all of these parameters are visible.

### PROCESS PARAMETERS

These parameters are used to control measuring process



#### LED measuring time

Controls the LED measuring time in older ALA instruments.

#### Measurement interval

Interval between the start of one measurement and the start of the next measurement. This is only active in the continuous mode. If the measurement process is longer than the interval, a continuous measurement is performed.

#### Number of chlorophyll (f) measurements

Number of measurements that are executed and averaged for one data set.

#### Measurement duration (f, fm, fo, transmission)

Please enlarge the measurement duration to enhance the accuracy of the measurements.

#### Minimum chlorophyll (f) concentration

Below this concentration no algae activity is calculated.

### Enter comment at start of measurement

If enabled, a comment related to the measurement can be entered. Please see window -> comment field for detailed configuration.

### Activity

Enables / disables the activity measurement (if available in the instrument)

### Continuous measurements

If enabled a continuous measurement is performed.

### Number of available sample channels

Shows the number of available sample channels – can only be changed by the bbe service

### Pump mode

Please select the way the sample pump is working.


ON: Always working

OFF: Off

Interval: Pump is working at the beginning of each measurement

### Pump time start up

Time, the pump is running before the first measurement. This is to get fresh sample and to remove remaining water from the measuring cell.



### Adaptation time (f, fo )


Period of time before the measurement starts. During this time, the LEDs are switched on to adapt the algae to the light conditions.


### Integration time fm measurement

Averaging time to smooth the fm results.

## ALARM PARAMETERS

These parameters control the alarm generation. Alarms are generated as soon as a threshold is exceeded. For each algae class and the total chlorophyll concentration a threshold can be set. In combination with a relais output, the alarms can be used to trigger an external device (like a sampler).






### Alarm limit

Set the threshold for the algae class

### Enable alarm

Enables the alarm.

### COMMON PARAMETERS





### Warm-up time

Period of time, before the measurement starts.

### Serial number


Serial number of the ALA.

### Software version

Shows the version of the firmware.

### FIT PARAMETERS

These parameters control the calculation of the results from the raw data.




### Name of class

Enter the name of the calibrated algae class or the name of the yellow substance.

### Enabled for fit

The enabled for fit section allows selecting maximum 5 of 8 calibrated algae classes for the calculation of the concentration of the different classes. Please take care, that all algae classes available in the water are activated.




### Yellow Substances Correction

Please select, whether this is a Yellow Substances fingerprint or not. Yellow substances content is not included in the total chlorophyll concentration.


### Offsets

The results of the offset calibrations are shown here.




## Transmission offsets / gradients

Shows the result of the transmission offset calibration.



## Fingerprints


Displays the results of the algae class calibrations.





## Standard deviations

Standard deviations of the fingerprints of each algae class. This standard deviation refers to variations within one specific algae class. The greater the standard deviation, the less important the measured value for this algae class



## Date of calibration

Shows the date of the last factory calibration.

## Global correction factor

Factor to adapt all measured chlorophyll results. All results will be multiplied with this factor.


## Cell factors

Factors to convert the chlorophyll concentration [ $\mu\text{g/l}$ ] results to cell counts [cell/l]. These factors can be defined separately for each algae class.

## F0 measurement brightness

Brightness of the LED light during the fo measurement in relation to the brightness during the standard chlorophyll measurement.



MEASUREMENT PARAMETERSDA values of the LEDs


Values result of the brightness regulation (service, do not change)

Temperature and Humidity offsets / gradients

Shows the result of the Temperature and humidity sensor calibration.

Required values of the LEDs

Values for the basic calibration of the ALA



### Temperature correction

Temperature correction factors for the given temperature ranges.

---

## OPERATING PRINCIPLES

---

### PRINCIPLES OF OPERATION

---

Due to the fact that algae of the same division contain a similar quantity and quality of photosynthetic pigments, their fluorescence excitation spectrum (with a fixed emission wavelength at 680nm) is significant. Thus, it is possible to differentiate divisions of algae by their fluorescence excitation spectrum. In addition to this, other fluorescing matter (for example, yellow substances) is detected to enhance the accuracy.

The bbe fluorometers use 6 LEDs for fluorescence excitation for algae differentiation. The LEDs emit light at 6 selected wavelengths (370nm, 470nm, 525nm, 570nm, 590nm and 610nm).

### DETERMINATION OF DIFFERENT ALGAE

The division of *chlorophyceae* (green algae) shows a broad maximum of fluorescence at the 470nm LED, which is caused by chlorophyll-a and -b. The *cyanophyceae* (blue-green algae) have their maximum at 610nm due to the photosynthetic antenna pigment *phycocyanin*. *Cyanophyceae* also contain chlorophyll-a if there is low intensity at 470nm. This is due to the masking effect of the *phycocyanin*. Furthermore, the high peak at the 525nm region for the *bacillariophyceae* originates from *xanthophyll fucoxanthin* and for the *dinophyceae* from *peridin*. The maxima at 470nm are caused by chlorophyll-a and -c. In our last analysed group, *cryptophyceae*, a significant maximum can be found at 570nm, which originates from *phycoerythrin*.

The different divisions of algae are first measured separately to calibrate the instrument. The measured spectra, or fingerprints, are then stored in the FluoroProbe.

During the measurement, the spectrum of the sample is loaded into the storage device of the instrument or sent to an external computer. The computer calculates the content of the different divisions of algae in the sample from the sample spectrum and the spectra of the separately measured algae divisions.

The concentration of every algae division is given in µg chlorophyll-a/l.

### DETERMINATION OF YELLOW SUBSTANCES

Yellow substances may fluoresce.

The UV LED (370nm) is used to measure yellow substances in the water. At 370nm it is possible to differentiate between algae (low signal) and yellow substances (high signal). The result of this measurement is given in relative units and not in weight/volume, because very different substances are detected. The result is mainly used to obtain a more accurate determination of algae classes – but it is also possible to determine the variations of the yellow substances.

### DATA ANALYSIS

#### Offset

The LEDs in the bbe fluorometer are switched on one after the other at high frequency. The fluorescence signal for each LED is taken and averaged during a given measuring time. The fluorescence values for each of the LEDs are given at the end. The concentration of the algae has to be calculated from these values.

The first step is to subtract the signal that does not come from the algae. This type of signal is also called offset. There are two types of offsets:

1. the signal from the bbe fluorometer itself, from the electronics and the optical system
2. the signal from fluorescent substances in the water that are not algae.

The first signal is determined by a measurement with distilled water – water without any fluorescent substance. The second signal can be determined by a measurement in ultra-filtrated sample water. This water only contains fluorescent substances that are not algae. Ultra-filtration is here understood as filtration with a 'mesh size' of 0.2 µm.

In the case of a varying concentration of these substances in the water, the subtraction fails. This is why the determination of yellow substances has been integrated into the bbe fluorometer. Yellow substances are treated as one algae class. Varying concentrations can be detected – assuming the fingerprint of the calibrated yellow substances matches the yellow substances in the sample.

The first signal only depends on the instrument and only varies within a relatively small range. So there is no need to calibrate these settings. The second signal depends on the sample, so there may be a need to calibrate this.

### Algae classes

After subtracting the offset from the measured fluorescence signals, the remaining signals have to be assigned to the different algae classes. A statistical calculation procedure in the software finds the best combination of concentrations of algae classes for the measured pattern of signals.

The total chlorophyll concentration is the sum of all detected concentrations of algae classes.

The signal can only be assigned to algae classes that are calibrated and activated in the instrument. Algae classes in the sample that are not part of the calibration will lead to miscalculations.

The different algae classes used for the statistical fit need to have significantly different fluorescence spectra such as the spectra calibrated by bbe in the bbe fluorometer. In case of queries, please ask **bbe Service**.

### MEASUREMENT OF THE ALGAE ACTIVITY

The fluorescence response to a very bright light impulse is measured under two different conditions. There is one measurement with and one without a very bright additional background light.

For the algae the light is the energy source and therefore like food. Without the additional background light, a large amount of the light impulse is used as "food". So hardly any light is emitted in the form of fluorescence. If there is enough light because of the background light, the level of the fluorescence response will increase considerably. If the cell is damaged, the cell will not use the light and even without additional background light there is a greater fluorescence response to the light impulse.

The so-called Genty parameter sums up the process and provides a measure of algae activity:

$$Genty = 100 * \frac{fm - fo}{fm} [\%]$$

This is:

Genty: shows the activity of the algae in percent whereby up to 75% is reached depending on the algae class and the physiological condition.

fo: fluorescence response without background light

fm: fluorescence response with background light

## MEASUREMENT OF FURTHER VALUES AND THEIR CALCULATION

### IN GENERAL

Generally, all measuring results are read into the PC as a figure between 0 and 4095. These have to be converted into the corresponding measuring values – for example, the transmission. As a unit of measurement for the figures, "digits" is used in the following.

The conversion of the figures into measuring values will generally be calculated by use of the following equation:

$$\text{Measurement Value} = \text{Gradient} * (\text{Figure} - \text{Offset})$$

The "offset" is the number which is read in if the measuring value is 0. "Gradient" is the conversion factor from the figures that are read into the measurement value.

Example:

0 digits is the same as  $-20^{\circ}\text{C}$

1200 digits is the same as  $100^{\circ}\text{C}$

So the offset is 200 [digits] and the gradient  $0.1 [^{\circ}\text{C}/\text{digit}]$ .

### TRANSMISSION

The transmission of the sample is measured at 5 different wavelengths: 470 nm; 525 nm; 570 nm; 590 nm; 610 nm.

The sensor is calibrated in a way that clear water has a transmission of 100% and "perfectly black water" has one of 0 %. The transmission can – in a limited way – also be used to compensate the fluorescence signal in cloudy water.

#### Calculation of the transmission from the sensor signals

$$\text{Transmission [\%]} = \text{Gradient [\%/digit]} * (\text{Signal [digit]} - \text{Offset})$$

#### Compensation of the chlorophyll concentration by the transmission

The compensation is based on measurements using clay as a model for turbidity. The coefficients are the same for all ALAs.

Two effects are compensated:

1. The higher offset due to additional reflection on the particles
2. The change of the signal (raw value) due to the mitigation of the light.

The correction is done separately for each LED.

Offset compensation:

$$\text{Offset corr. } (\lambda) = \text{Offset}(\lambda) + a * (100 - \text{Transmission } (\lambda))^2 + b * (100 - \text{Transmission } (\lambda)) * \text{Average Fingerprint (LED)}$$

a, b are parameters, determined by bbe.

Second step is the complete compensation of the raw value:

$$\text{Raw corr. } (\lambda) = \text{Raw } (\lambda) \cdot (1 + \text{Factor}_{\text{mitigation}}(\lambda) \cdot (100 - \text{Transmission}))$$

After doing the correction, the calculation of the chlorophyll concentrations is done as normal.

### TEMPERATURE OF THE SENSOR

The very small temperature drift of the sensor can be compensated additionally. Compensation factors for pre-set temperature ranges can be set.

#### Calculation of the temperature compensation

$$\text{Chlorophyll concentration } [\mu\text{g/l}] = \text{measured chlorophyll concentration } [\mu\text{g/l}] \cdot \text{correction factor (Temp)}$$

If the measured chlorophyll concentration at an internal temperature of 20°C is 100µg/l and at an internal temperature of 30°C is 101µg/l, the correction factor is 0.99 at 30°C.

### BRIGHTNESS OF THE LED

The fluorescence signal depends directly on the brightness of the excitation lamp of the diodes. The brightness is regulated to the target value at the beginning of each measurement. Afterwards the brightness and the correction values are determined regularly.

## MEASURING PROCEDURE

The bbe Algae Online Analyser was calibrated before delivery with bbe standard algae types and yellow substances. In case of the presence of algae types that are not calibrated or in case of different types of yellow substances, the bbe Algae Online Analyser needs to be calibrated to these special needs first. This is most important if there are only low concentrations of algae (< 10µg/l).

### INITIAL MEASURING PROCEDURE

- Start the instrument and in case of an external PC also the PC and the AOA Software
- Connect the sample tubes to the sample supply and test it by choosing Test - Pump
- Start the measurement (Measurement – Start ).
- A window appears where the COM port has to be chosen.
- At the end of the measurement, stop the measurement. The data are saved automatically in the database .

### OFFSETS

The following possibilities are available to deal with the offset:

1. The offsets are most important if the concentration of chlorophyll is low, for example lower than 10µg/l. In such cases, the fingerprint of the yellow substances should be recalibrated by using ultra-filtrated sample water (see below).
2. In case of stable yellow substance concentrations, the results may become more stable if the yellow substances are deactivated. In this case, the offset “ultra-filtrated water” is used for the calculation. Please calibrate this offset before choosing this option.

## MEASURING PROCEDURE FOR LOW CHLOROPHYLL CONCENTRATION

The lower the chlorophyll concentration, the more exact the offsets and the yellow substance offset have to be.

## BACKGROUND FLUORESCENCE

The first consideration is whether the background fluorescence (fluorescence of the ultra-filtrated sample) varies its concentration or not.

- *In case of a stable background fluorescence, the best choice is to deactivate the yellow substances in the "Parameter of fit" window and perform an offset calibration with ultra-filtrated sample water.*
- *In case of varying background fluorescence, the best choice is to measure this fluorescence as yellow substances. Please perform a fingerprint calibration of the yellow substances with ultra-filtrated sample water. The yellow substances need to stay activated.*

## CLEANING

Please ensure that all tubes and the measuring chamber are really clean. Especially after measuring high concentrations of chlorophyll, multiple rinsing is recommended. To check the cleaning status of the whole system, please perform a measurement with distilled water. The result should be 0. If not, take also in consideration that the tubing might be dirty inside.

Please note: Deionized water may also contain fluorescing substances.

## MEASURING TIME

For low concentrations of chlorophyll, please prolong the measuring time to 60s for lower noise.

## SEDIMENTATION

In case of sedimentation of algae in the sensor, choose the operation mode "Always on" for the internal pump.

## SAMPLE PREPARATION

To dilute a sample, please use algae nutrient solution. Tap water may lead to a different status of the algae. The prepared sample is also more stable in a nutrient solution.

Please note: algae may change their fluorescence properties in relatively short time periods.

In case of deactivated yellow substance fingerprints, the dilution has to be carried out with ultra-filtrated sample water.

## CALIBRATION PROCEDURE

---

### GENERAL CONSIDERATIONS:

- To add any type of solution to the chlorophyll sensor, please use the syringe supplied as an accessory. Disconnect the internal pump and connect the syringe at the sample outlet. To apply a solution, suck it up into the chlorophyll sensor by using the syringe.
- Make sure that the chlorophyll sensor is really clean! Do a cleaning procedure first.

### CALIBRATION OF THE OFFSETS:

- For the calibration of the offsets, please use ultra-filtrated sample water and distilled water.
- Apply distilled water. Start using Calibration -> Offsets (distilled).
- Apply the values after the calibration.
- Apply ultra-filtrated sample water and start using Calibration -> Offsets (ultra-filtrated) .
- Apply the values after the calibration.

### CALIBRATION OF THE ALGAE FINGERPRINTS:

- Do an offset calibration first.

- For the calibration of the algae fingerprints please use pure solutions of the algae division with a known chlorophyll concentration. The concentration should be in the range of 50-100µg/l.
- Apply the algae solution.
- Choose -> Calibration -> Fingerprints in the software, choose the algae class and enter the concentration.
- Choose the type of dilution water used for the calibration.
- Start the calibration.
- Apply the values after the calibration.

#### CALIBRATION OF THE YELLOW SUBSTANCES:

- Do an offset calibration first.
- For the calibration of the yellow substances please use ultra-filtrated sample water.
- Apply the solution.
- Choose Calibration -> Fingerprints in the software, choose the yellow substances and enter 1 r.u. as the concentration and start the calibration.
- Apply the values after the calibration.

#### CALIBRATION BY FACTORS:

The assumption for this type of calibration is that there is a number of results of the Algae Online Analyser and of a reference method that show a certain factor for one algae class. It is highly recommended to use only measurements where this algae class is dominant and the concentration at least 10µg/l.

- Choose -> Calibration -> Correction factor in the software
- Choose the algae class and enter the determined factor.
- Apply the values.

#### CALIBRATION OF FINGERPRINTS WITH SUBSEQUENT WET-CHEMICAL ANALYSIS:

If the Algae Online Analyser is calibrated with a reference solution where the concentration is determined subsequently, the following steps have to be carried out:

- Do a fingerprint calibration as described above, but enter 1 for the concentration of the reference solution
- Do a wet-chemical analysis
- Use the "Calibration by factor" feature to adapt the calibration to the result of the subsequent wet-chemical analysis.

Attention: the Algae Online Analyser is not ready for use until the last step has not been carried out.

#### BATCH MODE OPERATION


---

The batch mode can be used to operate the instrument with single samples for example in the lab. While in the online mode everything is automatically done by the software, in the batch mode every step can be influenced by the user.

To switch in the batch mode, activate "batch" in the view menu. The batch panel to operate the instrument is displayed. Depending on the features of the chlorophyll sensor the Activity and Transmission click boxes are greyed or not greyed.

The sample can be taken from a canister or similar container. To get the sample into the sensor, either the internal pump or a syringe can be used.





Single channel version



Multi-channel version

To perform a measurement, please carry out the following steps:

1. Mount the tubes on the instrument so that the container and effluent can easily be reached.
2. Put the suction tube in the container and click "Pump". The pump runs as long as the value given in "Fill sensor" in the Parameters. Hint: The pump should work as long as it takes to pump sample out of the effluent of the sensor.

In case of a multi-channel instrument 2 or more pump buttons are available. Please use the button of the appropriated sample.

3. Press "Start" to start the measurement. If Transmission measurement and/or cleaning shall be executed at the end of the measurement procedure, click the according checkboxes. The measurement times etc. are given in the Parameters menu.
4. Put the suction tube in clean water and press pump to rinse the sensor. If there are variations in the different samples, this has to be carried out very carefully, i.e. 3-4 times.

It also helps to use the clean option while pure water is in the chamber. This requires one more rinsing afterwards. If there is enough sample water, it is recommended to use the next sample for rinsing.

## DATA EVALUATION

---

---

To evaluate the data, either the PC of the instrument itself can be used or any other PC on which the AOA software is installed.

On a desktop PC, please install the bbe++. Now any bbe++ database file (\*.BDB) can be opened for evaluation. The file can either be on the PC or in the LAN network.

## SERIAL DATA EXCHANGE

---

To change settings for the serial data exchange go to Tools -> settings -> Text results.

The protocol is ASCII-based. The Algae Online Analyser waits for a character and then sends an answer if a command has been recognised. The following commands are allowed:

- "h" request for the header with the names of the single data.
  - "u" request for the units of the single data.
  - "d" request for the actual dataset.
  - "s" to start the measurement. In this case there will be no answer string.
  - "p" to stop the measurement. In this case there will be no answer string.
- The answer is a string terminated with CR and LF. Each string contains a fixed number of separated fields. The data fields and setting are defined in the output view chosen in Tools -> setting -> text results.

The commands „h“ and „u“ always return header and units. Please note the unit string contains the character “μ” ASCII code 230. This might appear as „æ“ when ANSI coding is applied.

The “d” command always returns the actual dataset. If there is no new dataset since the last request, the same dataset is delivered.

The first 5 data fields are fixed and always transmitted:

Feld	Typ	Inhalt
1	Text	Checksum in the format "Pxx", where xx is a hexadecimal value in text representation for example "P0F" - P: marker for checksum; 0F: hexadecimal representation of 15.
2	Text	Number of software version
3	Zahl	Number of data fields
4	Text	Day, month, year in the format "DD.MM.YYYY"
5	Text	Date and time in the format "hh:mm:ss"

The first field is the checksum of the whole string. The format of the checksum is a P (capital) and a two-digit hexadecimal number. The checksum is an XOR function of the following characters, excluding the TAB directly after the checksum and excluding the CR and LF. All the other TABs are included. The calculation of the checksum starts with 0. The 0 is linked by XOR to the ASCII number of the first character after the first TAB. The result is linked to the ASCII number of the second character, and so on.

“C” example for the calculation of the checksum:

```
char * response; /* response enthält den vollständigen Antwortstring */
int i = 4; /* Prüfsumme und erstes TAB nicht berücksichtigen */
unsigned char checksum = 0; /* Startwert ist 0 */

while (response[i] && (response[i] != '\r') && (response[i] != '\n'))
{
    checksum = checksum ^ response[i];
}
```

```

    i = i + 1;
}

```

The value in the checksum has to be the same as the two-digit number after the P.

Examples: The data fields might vary according to the settings:

External device send the request "h". Answer:

```

P1B  02.00 33   Datum Uhrzeit   Konz. 1   Konz. 2   Konz. 3   Konz. 4
      Konz. 5   Konz. 6   Konz. 7   Konz. 8   Gesamtkonz. Transmission
      Tiefe Temp. Probe      Zellz. 1   Zellz. 2   Zellz. 3   Zellz. 4
      Zellz. 5   Zellz. 6   Zellz. 7   Zellz. 8   Gesamtzellzahl   Trans. 700
nm   LED 3 LED 4 LED 5 LED 6 LED 7       LED 8 Druck Temp. LEDs   Temp. Sensor

```

Externer Rechner sendet „u“, Antwort:

```

P4D  02.00 33   Datum Uhrzeit   µg/l µg/l µg/l µg/l µg/l µg/l µg/l
      µg/l %     m °C   Zellen/ml Zellen/ml Zellen/ml Zellen/ml
      Zellen/ml Zellen/ml Zellen/ml Zellen/ml Zellen/ml Digits
      Digits   Digits   Digits   Digits   Digits   Digits   bar
      °C      °C

```

Externer Rechner sendet „d“, Antwort:

```

P33  02.00 33   25.11.2014  14:27:38   0,00 1,66 0,00 0,00 0,00 0,00 0,00
      0,00 1,66 100,00   -0,10 15,10 0 1656 0 0 0 0
      0 1656 354,22   4,29 4,01 12,29 2,37 7,68 2,59 0,99 16,95 15,46

```

## TECHNICAL DATA

Description	Value
Measurands	total chlorophyll [ $\mu\text{g Chl-a/l}$ ]
	concentration of green algae [ $\mu\text{g Chl-a/l}$ ]
	concentration of cyanobacteria [ $\mu\text{g Chl-a/l}$ ]
	concentration of diatoms [ $\mu\text{g Chl-a/l}$ ]
	concentration of cryptophytes [ $\mu\text{g Chl-a/l}$ ]
	yellow substances
	photosynthetic activity (Genty) – Option
	transmission (at 5 wavelengths)
	water temperature
Chlorophyll	0 - 200 $\mu\text{g Chl-a/l}$
Measurement procedure	spectral fluorometry
Resolution	0.01 $\mu\text{g Chl-a/l}$
Transmission	0 - 100 %
Measuring chamber cleaning	automatic cleaning piston
Housing material	V4A steel   aluminium   coated steel plate
Weight	19 kg
Dimensions (H x W x D)	420 x 600 x 200 mm
Protection class	IP 54
Mains voltage	110/ 240V ; 50/60Hz
Power consumption	100 W
Fuse	2 x 3,15A T
Sample temperature	0 – 40° C
Sample volume	30 ml
Temperature range while on transport	0 – 50° C
Relative humidity	Up to 95%, not condensing
Sunlight	indirect
Maintenance interval	> 7 days
Sample inflow	free flow / tube pump
Connection Inlet	6 mm internal hose connection
Connection outlet	6 mm internal hose connection
Sample pump	Marprene Tubing ID = 4.8 mm
Valves (option)	Silicon tubing ID=4.80 / OD=7.90mm
PC	ATOM N270 12" touchscreen panel - ca. 1.6 GHz, 1 GB RAM
OS	Windows XP Pro (SP3)

---

Description	Value
Hard disk	80 GB
Outputs (options)	Modem
	2 x analogue output 4-20 mA
	2 x relay output
	SDI-12 with bbe converter

## MAINTENANCE

---

### WEEKLY MAINTENANCE :

---

#### MAINTENANCE WORK THAT HAS TO BE DONE:

- ⇒ Check the inflows and drains for blockages
- ⇒ Check the instrument for leaks
- ⇒ Clean the hose

### MONTHLY MAINTENANCE

---

#### MAINTENANCE WORK THAT HAS TO BE DONE:

- ⇒ Exchange the pump hoses
- ⇒ Clean the measuring cuvette


### HINTS FOR MAINTENANCE WORK

---

#### CLEANING OF THE HOSE SYSTEM

1. Stop the measurement
2. Open the hose saddle of the sample pump
3. Connect the water-filled syringe to the instrument instead of the sample
4. Rinse the hose system with pump using the syringe

#### EXCHANGE THE PUMP HOSES



The pump turns clockwise. Make sure that the outlet at the right side of the pump is connected to the inlet of the sensor at the bottom of the sensor!

#### Hose diameter and hose types


In any case, only hoses in accordance with the specifications must be used!

Marprene hose 4.8 x 1.6 mm (inner diameter x wall thickness)

The inner diameter can be varied according to the flow needed. To vary the diameter, please refer to the chapter 'pump'.

### CLEAN THE MEASURING CUVETTE

1. Stop the measurement
2. Switch the instrument off
3. Disconnect influent and effluent tubes directly at the sensor
4. Disconnect the electrical connector of the sensor
5. Open the mounting screws at the top and the bottom of the sensor
6. Take out the sensor and turn it upside down
7. Open the 4 small screws around the inlet of the sensor
8. Take off the lid
9. Clean the cuvette inside and the lid with a smooth cloth
10. Mount the lid again and be careful with the seal
11. Mount everything again vice versa



View from the bottom of the sensor. Arrows mark the socket screws which have to be removed to demount the removable measuring chamber base (lower lid).



## CHANGE FUSES

---

If the instrument does not start, a fuse may have blown. There are two fuses on the right side of the instrument. To change them:

- Disconnect the mains
- Open the lid
- Take out the fuse
- Replace it (3.1 A T)
- Close the lid

---

## FAULT FINDING

---

---

### THE ALGAE ONLINE ANALYSER WILL NOT START

---

- Is the Algae Online Analyser connected to the mains?
- Is the Algae Online Analyser switched on?
- Are the fuses OK?

---

### PUMP IS NOT RUNNING

---

- Test the pump by using the “Test pump” option.
- Is the pump mode set to “Always off”

---

### NO ALGAE CONCENTRATION SHOWN

---

- There is no sample in the sensor. Maybe because inlet and outlet of the sensor are mismatched
- No sample supply
- Is the pump mode set to “Always off”

---

### GENTY RESULTS ARE ZERO

---

- The concentration is lower than the minimum concentration

---

### RESULTS TOO LOW

---

- Measuring cuvette is dirty
- The cleaning device is not working. Please check

---

### RESULTS TOO LOW OR TOO HIGH

---


- You are using a different method of chlorophyll calibration than HPLC. So a certain factor has to be used

---

### HUMIDITY DETECTED IN SENSOR

---

- The software shows the following message:



Please disconnect the sensor from the power supply and the sample and contact bbe service.

## TROUBLESHOOTING

---

---

Please call us:

Tel.: +49 431 380 40-40 (Service)

Tel.: +49 431 380 40-0

Fax: +49 431 380 40-10

Or send an email to:

**e-mail:** [bbe@bbe-moldaenke.de](mailto:bbe@bbe-moldaenke.de)

[www.bbe-moldaenke.de](http://www.bbe-moldaenke.de)

We will gladly help you.